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INSTRUCTION FOR USE

1 IDENTIFICATION OF THE DEVICES

These Instructions for Use (IFU), provided in hard copy, accompany Class Ila surgical instruments manufactured by BRM Extremities S.r.l. and intended for implantation of BRM Extremities Bioplan device. These devices, listed in Table 1, are considered accessories to the BRM Extremities Bioplan implantable device and may be used only for activities related to the implantation and explantation of the indicated device.

TYPE OF INSTRUMENT	CODE	COMMERCIAL NAME
Trial	SBP110P007	Bioplan - Trial Size 7
	SBP110P008	Bioplan - Trial Size 8
	SBP110P009	Bioplan - Trial Size 9
	SBP110P010	Bioplan - Trial Size 10
	SBP110P011	Bioplan - Trial Size 11
	SBP110P012	Bioplan - Trial Size 12
	SBP110P013	Bioplan - Trial Size 13

Table 1 – Class Ila surgical instruments accessories to the Bioplan device.

For any clarifications or requests, please contact the manufacturer at the following numbers:

- BRM Extremities S.r.l.**
- Registered Office: Via Lorenzo Mascheroni, 29 - 20145 Milano (MI) - Italy
 - Operational Headquarters: Via Privata ai Campi, 7 - 23862 Civate (LC) – Italy
 - info@brm-extremities.com
 - Tel. +39 0341 1693087

1.1 GENERAL INFORMATION

These IFU refer to Trials, surgical instruments that cannot be connected to active devices, intended to be used for the correct implantation of Bioplan BRM Extremities devices. The devices are supplied in NON-STERILE form and are reusable, subject to conditioning. Please refer to the implantable device IFU and surgical technique for a list of all instruments to be used for implanting the device.

2 PRELIMINARY WARNINGS

Carefully read and follow the warnings below. Failure to comply with the warnings described in these IFUs will result in the immediate invalidation of any warranty on the medical device. BRM Extremities Srl is not responsible for any damage to people or things resulting from failure to comply with the warnings.

BEFORE USE, READ AND UNDERSTAND THE INSTRUCTIONS FOR USE (IFU) ENTIRELY. FAILURE OR NEGLIGENCE IN COMPLYING WITH THE INSTRUCTIONS FOR USE COULD CAUSE THE DEVICE TO MALFUNCTION, WITH POSSIBLE DAMAGE AND INJURY TO THE USER.

DEVICES MUST BE CHECKED, WASHED AND STERILIZED BEFORE USE, AS DESCRIBED IN SECTION 6.

FOLLOW THE INSTRUCTIONS FOR USE DESCRIBED IN THE IFU. USES DIFFERENT THAN THOSE INTENDED COULD EXPOSE THE PATIENT OR USER TO UNCONSIDERED DANGERS.

DO NOT MODIFY HE DEVICE IN ANY WAY.

THE DEVICES ARE INTENDED FOR USE EXCLUSIVELY BY ORTHOPAEDIC SURGEONS WITH CLINICAL EXPERIENCE IN THE IMPLANT SURGICAL TECHNIQUE.

FOLLOW THE SURGICAL TECHNIQUE INDICATED BY THE MANUFACTURER AND USE INSTRUMENTS ONLY FOR THEIR INTENDED PURPOSE.

CHECK THE INTEGRITY OF THE DEVICE. IN CASE OF DAMAGE, DO NOT USE AND CONTACT THE MANUFACTURER FOR REPLACEMENT.

IT IS STRICTLY FORBIDDEN TO USE THESE INSTRUMENTS IN COMBINATION WITH OTHER DEVICES/ACCESSORIES NOT MENTIONED IN THESE INSTRUCTIONS (SECTION 8), AS GEOMETRIC-DIMENSIONAL COMPATIBILITY MAY NOT BE ENSURED.

3 INTENDED USE

The intended use of the devices is the implantation of BRM Extremities Bioplan devices, according to the surgical technique defined by the manufacturer and in conjunction with the other components of the instrumentation.

REFER TO THE SURGICAL TECHNIQUE FOR THE IMPLANTATION OF THE BIOPLAN DEVICE IN ORDER TO IDENTIFY THE STEP OF USE OF THESE DEVICES.

3.1 INTENDED PATIENT GROUP

These devices are accessories of Bioplan implantable device. For identification of the intended patient group, please refer to the IFU of this implantable device.

3.2 INDICATION AND EXPECTED CLINICAL BENEFITS

The surgical instruments in question are accessories to the implantation of the Bioplan device. It is not possible to define a clinical benefit for an instrument taken individually. These instruments are part, together with others belonging to classes I and/or II, of a set whose clinical benefit consists in the simplification of the related implantation and/or explantation operations. The expected clinical benefit can therefore only be assessed in connection with the implantation of the medical device.

3.3 CONTRAINDICATION

These devices are accessories of the Bioplan implantable device. For identification of contraindications, please refer to the IFU of the implantable device.

3.4 RESIDUAL RISKS AND SIDE EFFECTS

Given the intended use and the materials used to manufacture the devices, the development of allergic reactions related to the chemical elements in the metal alloys is extremely rare. In the case of particularly sensitive individuals, contact the manufacturer. No other side effects related to these devices are known.

Intra-operative complications associated with the surgical procedure as a whole, of which the individual instruments in this IFU are an integral part, are reported in the IFU of the individual implantable device.

The following residual risks related to individual instruments are also identified:

- Failure to apply the manufacturer's intended surgical technique could result in injury to the patient, including permanent injury.
- In case of complete loss of function of an instrument during the performance of surgery, e.g., due to a fall, it is the surgeon's responsibility to consider proceeding with alternative modalities or discontinuing the surgery.
- The indicated reprocessing method may not be suitable for complete prion inactivation. Follow local or national guidelines.

4 INFORMATION ON IMPLANTABLE DEVICES

N.A. as the devices in question are not implantable devices.

5 DEVICE DESCRIPTION AND EXPECTED PERFORMANCE

5.1 TRIALS

Trials are surgical instruments for manual use with 'AO' type handles. They have a cannulation that allows them to be inserted on a guide pin. They consist of a shank, on which reference notches are marked to determine their positioning, and an end that is geometrically and dimensionally similar to that of the corresponding implantable device. The trials are inserted into the anatomical area to select the size of the Bioplan device that best fits the patient's anatomical features as well as to determine its positioning.

They are available in one variant, in seven sizes, equal to the available Bioplan device sizes.

5.2 MATERIALS CONSTITUTING THE DEVICE

Trials are made of AISI 630 stainless steel.

5.3 PACKAGING AND STERILITY

The devices are supplied non-sterile and individually packaged in transparent polyethylene bags.

TAKE THE DEVICE OUT OF THE PACKAGING BEFORE PROCEEDING WITH REPROCESSING OPERATIONS

IF THE PACKAGING IS DAMAGED, DO NOT USE THE DEVICE AND CONTACT THE MANUFACTURER

6 DEVICE CONDITIONING

Devices should be cleaned and sterilised before use, following instructions in the document: **Instructions for Care and Conditioning of Medical Devices and Accessories intended to be Conditioned by the User** attached to this IFU. If the conditioning instructions are lost, contact the manufacturer.

7.1 PREPARING THE DEVICE

After being sterilized, the devices can be used for the application of interest. Their handling is reserved to qualified personnel.

EXERCISE EXTREME CAUTION TO AVOID CONTAMINATION OF THE DEVICE PRIOR TO USE. IF THE DEVICE IS SHOWN OR A LOSS OF STERILITY IS POSSIBLE, THE DEVICE MUST BE REPROCESSSED. THE USE OF POTENTIALLY CONTAMINATED DEVICES COULD CAUSE INFECTION, RESULTING IN THE NEED FOR MEDICAL TREATMENT AND POSSIBLE IMPLANT FAILURE.

7.2 MODE OF USE

The devices are intended for the implantation of Bioplan device, according to the relevant surgical technique. The operating methods of the individual instruments are summarized below.

BEFORE USE, CHECK THAT THE DEVICES SHOW NO SIGNS OF DAMAGE OR CORROSION AND THAT THE MARKINGS ARE CORRECTLY LEGIBLE. IF THIS IS NOT THE CASE, CONTACT THE MANUFACTURER TO REPLACE THE DEVICES.

7.2.1 TRIALS

- Take the trial size consistent with the size of the device to be implanted.
- Connect the trial to the handle.
- Insert the trial into the implant area and check if the chosen size is appropriate. Remove the trial and disconnect it from the handle.
- If necessary, repeat the operation with trials of other sizes, to identify the most suitable for the patient.

BEFORE USE, CHECK THAT THE TRIAL IS CORRECTLY ATTACHED TO THE HANDLE.

8 DEVICES INTENDED TO BE USED TOGETHER

The instruments to be used together are summarized in the following table.

DEVICE	DEVICES INTENDED TO BE USED TOGETHER
SBP110P007	Bioplan - Trial Size 7
SBP110P008	Bioplan - Trial Size 8
SBP110P009	Bioplan - Trial Size 9
SBP110P010	Bioplan - Trial Size 10
SBP110P011	Bioplan - Trial Size 11
SBP110P012	Bioplan - Trial Size 12
SBP110P013	Bioplan - Trial Size 13

9 TREATMENT AFTER USE

After use, devices should undergo appropriate reconditioning and inspection activities as described in the **Instructions for Care and Conditioning of Medical Devices and Accessories intended to be Conditioned by the User**.

10 SHELF LIFE

The devices are supplied as non-sterile versions and do not have a defined expiration date. At each reprocessing, however, it is necessary to verify that the devices:

- have no visible signs of wear, breakage or deformation;
- that the connection with the handle or motor is stable, in good condition and working properly.

Precise directions regarding these activities can be found in the **Instructions for Care and Conditioning of Medical Devices and Accessories intended to be Conditioned by the User**.

11 DISPOSAL OF THE DEVICE AND PACKAGING

To dispose of the device, follow the waste management regulations of the Local Management Authority (Special hospital waste). Do not dispose of the product or its packaging in the environment under any circumstances.

Packaging	Material family	Decision 97/129/CE
Transparent bag	Polyethylene	C/DP81
IFU	Paper	PAP22

12 TRANSPORT AND STORAGE CONDITIONS

The devices can be transported and stored at room temperature and must be kept out of direct light.

13 INCIDENT REPORTING AND SURVEILLANCE ACTIVITIES

ATTENTION: As per current legislation, the user is required to report any serious incident occurring in relation to the device to the manufacturer and the competent authority of his country.

It is also good practice to inform the manufacturer of any minor incidents, complaints, suggestions, simple feedback or possible improvements to the device. The information collected will be evaluated within the Surveillance Plan implemented by the manufacturer and will contribute to improving the product.

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TIPOLOGIE	CODE	NON COMMERCIAL
Tailles d'essai	SBP110P007	Bioplan - Trial Size 7
	SBP110P008	Bioplan - Trial Size 8
	SBP110P009	Bioplan - Trial Size 9
	SBP110P010	Bioplan - Trial Size 10
	SBP110P011	Bioplan - Trial Size 11
	SBP110P012	Bioplan - Trial Size 12
	SBP110P013	Bioplan - Trial Size 13

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