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Tabella 1 – Reusable surgical instruments for the implantation of Bioplan device.

For any clarifications or requests, please contact the manufacturer at the following numbers:

- BRM Extremities S.r.l.**
• Registered Office: Via Lorenzo Mascheroni, 29- 20145 Milano (MI) - Italy
• Operational Headquarters: Via Privata ai Campi, 7- 23862 Civate (LC) – Italia
• info@brm-extremities.com
• Tel. +39 0341 1693087

1.1 GENERAL INFORMATION

These IFUs refer to reusable surgical instruments intended to be used for the implantation and explantation of BRM Extremities Bioplan devices, that cannot be attached to active devices.

The devices are supplied in NON STERILE form and are reusable, subject to conditioning.

The devices can be divided into the following categories:

- Screwdriver
- Extractor Screwdriver
- Probe
- Guide wire Ø 2 mm

Refer to the implantable device IFUs and surgical technique for a list of all instruments to be used for anatomic area preparation and implantation of devices.

2 PRELIMINARY WARNINGS

Carefully read and follow the warnings given below. Failure to follow the warnings described in these IFUs will result in the immediate voiding of any warranty on the medical device. BRM Extremities S.r.l. shall not be liable for any damage to persons or property because of failure to follow the warnings.

BEFORE USE, READ AND UNDERSTAND THE OPERATING INSTRUCTIONS (IFU) IN THEIR ENTIRETY. FAILURE OR NEGLECT TO FOLLOW THE INSTRUCTIONS FOR USE COULD RESULT IN MALFUNCTION OF THE DEVICE, POSSIBLY CAUSING DAMAGE AND INJURY TO THE USER OR THE PATIENT.

DEVICES MUST BE CHECKED, WASHED AND STERILIZED BEFORE USE, AS DESCRIBED IN SECTION 6.

FOLLOW THE INDICATION DESCRIBED IN THE IFU. USES OTHER THAN THOSE INTENDED COULD EXPOSE THE PATIENT OR USER TO UNFORESEEN DANGERS.

DO NOT MODIFY THE DEVICE IN ANY WAY.

THE DEVICES ARE INTENDED FOR USE ONLY BY ORTHOPEDIC SURGEONS WITH CLINICAL EXPERIENCE ON SURGICAL IMPLANTATION TECHNIQUE.

FOLLOW THE SURGICAL TECHNIQUE INDICATED BY THE MANUFACTURER AND USE INSTRUMENTS ONLY FOR THEIR INTENDED PURPOSE.

CHECK THE INTEGRITY OF THE DEVICE. IF DAMAGED, DO NOT USE AND CONTACT THE MANUFACTURER FOR REPLACEMENT.

IT IS STRICTLY FORBIDDEN TO USE THESE INSTRUMENTS IN CONJUNCTION WITH OTHER DEVICES/ACCESSORIES NOT INDICATED IN THESE INSTRUCTIONS (SECTION 8). AS GEOMETRIC-DIMENSIONAL COMPATIBILITY MAY NOT BE ENSURED, IT IS FORBIDDEN TO USE THE DEVICES IN CONNECTION WITH ACTIVE DEVICES.

3 INTENDED USE

The intended use of the devices is the implantation and explantation of BRM Extremities Bioplan devices, according to the surgical technique defined by the manufacturer and in conjunction with the other components of the instrumentation.

REFER TO THE SURGICAL TECHNIQUE FOR IMPLANTATION OF THE BIOPLAN DEVICE IN ORDER TO IDENTIFY THE STAGE OF USE OF THESE DEVICES

3.1 INTENDED PATIENT GROUP

These devices are accessories of the implantable Bioplan device. Please refer to the IFUs of the reference implantable device for identification of the intended patient group.

3.2 INDICATION AND EXPECTED CLINICAL BENEFITS

The surgical instruments are accessories to the implantation of Bioplan device. It is not possible to define a clinical benefit for an instrument taken individually. These instruments are part, together with others, belonging to classes and/or to the, of a set whose clinical benefit consists in the simplification of the related implantation and/or explantation operations. The expected clinical benefit is therefore assessable only in relation to the implantation of the medical device.

3.3 CONTRAINDICATION

These devices are accessories of the implantable Bioplan device. Please refer to the IFUs of the reference implantable device for identification of the contraindications.

3.4 RESIDUAL RISKS AND SIDE EFFECTS

Given the intended use and the materials used to manufacture the devices, the development of allergic reactions related to the chemical elements in the alloys is extremely rare. In the case of particularly sensitive individuals, contact the manufacturer. No other side effects related to these devices are known.

Intra-operative complications associated with the surgical procedure as a whole, of which the individual instruments in this IFU are an integral part, are reported in the IFU of the individual implantable device.

The following residual risks related to individual instruments are also identified:

- Failure to apply the manufacturer's intended surgical technique could result in injury to the patient, including permanent injury.

- In case of complete loss of function of an instrument during the performance of surgery, for example due to a fall, it is the surgeon's responsibility to consider proceeding with alternative modalities or discontinuing the surgery.
- The indicated reconditioning method may not be suitable for complete prior inactivation. Follow local or national guidelines.

4 INFORMATION ON IMPLANTABLE DEVICES

N.A. as the devices in question are not implantable devices.

5 DEVICE DESCRIPTION AND EXPECTED PERFORMANCE

5.1 SCREWDRIVER

The insertion screwdriver is a cannulated steel tool with an 'AO' connection for connecting to the handle (see point 8 for details) and a hexagonal bit for coupling with the Bioplan device. It has numbered marks on the cylindrical body. It is used for inserting the Bioplan device into the implant site.

5.2 EXTRACTOR SCREWDRIVER

The extractor screwdriver is a steel tool equipped with an 'AO' connection for connection to the handle (see point 8 for details) and a tip that allows coupling with grooves on the Bioplan device. It is used to extract the Bioplan device from the implant site.

5.3 PROBE

The probe is a cannulated instrument with an 'AO' connection used for detecting the sinus of the tarsus.

5.4 GUIDE WIRE

The Guide Pin is a 2 millimetre diameter blunt-tipped wire used to provide guidance and support during implantation of the Bioplan device.

5.5 CONSTITUTENT MATERIALS OF THE DEVICE

The guide pin is made of AISI 316L stainless steel, the other devices of AISI 630.

5.6 PACKAGING AND STERILITY

The devices are supplied non-sterile and individually packaged in transparent polyethylene pouches.

REMOVE THE DEVICE FROM THE PACKAGING BEFORE PROCEEDING WITH REPROCESSING OPERATIONS

IF THE PACKAGING IS DAMAGED, DO NOT USE THE DEVICE AND CONTACT THE MANUFACTURER

6 DEVICE CONDITIONING

Devices should be cleaned and sterilized before use, following instructions in the document *Recommended processing instructions for medical devices and accessories intended to be conditioned by the user* attached to this IFU. If the conditioning instructions are lost, contact the manufacturer.

7 USE OF THE DEVICE

7.1 PREPARATION OF THE DEVICE

After being sterilized, the devices can be used for the application of interest. Their handling is reserved for qualified personnel.

OPERATE WITH EXTREME CAUTION TO AVOID CONTAMINATION OF THE DEVICE BEFORE USE. IF THE DEVICE IS DROPPED OR POSSIBLE LOSS OF STERILITY, THE DEVICE MUST BE REPROCESSED. IMPLANTATION OF POTENTIALLY CONTAMINATED DEVICES COULD CAUSE INFECTION, RESULTING IN THE NEED FOR MEDICAL TREATMENT AND POSSIBLE IMPLANTATION FAILURE.

7.2 MODE OF USE

The devices are intended for implantation of the Bioplan device, according to the relevant surgical technique. The modes of operation of individual instruments are summarized below.

BEFORE USE, CHECK THAT THE DEVICES SHOW NO SIGNS OF DAMAGE OR CORROSION AND THAT THE MARKINGS ARE PROPERLY LEGIBLE. IF NOT, CONTACT THE MANUFACTURER FOR REPLACEMENT OF THE DEVICES.

7.2.1 SCREWDRIVER

- Connect the screwdriver to the handle.
- Insert the screwdriver on the guide pin and bring the tip into contact with the head of the Bioplan device.
- When the tip of the screwdriver is firmly engaged in the screw head, start rotation clockwise.
- Position the Bioplan device, using the notches on the screwdriver body as a reference.
- Remove the screwdriver.

7.2.2 EXTRACTOR SCREWDRIVER

- Connect the screwdriver to the handle
- Bring the tip of the screwdriver into contact with the head of the Bioplan device, coupling it with the grooves on the body of the implanted device.
- When the screwdriver tip is firmly engaged, start rotation counter-clockwise and withdraw the device.

THE USE OF SCREWDRIVERS OTHER THAN THOSE PRESCRIBED BY THE MANUFACTURER COULD RESULT IN DAMAGE TO THE HEAD AND BODY OF THE BIOPLAN DEVICE, POTENTIALLY MAKING IT IMPOSSIBLE TO COMPLETE IMPLANTATION OR REMOVAL.

7.2.3 PROBE AND GUIDE WIRE

- Connect the probe to the handle.</