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- Handles: Manual tools equipped with a handle and an attachment to accommodate an instrument.

- Kirschner wires: Perforated cylindrical containers designed to contain the Kirschner wires.

- Spingifilo: Manual instrument that can be used with Silktote Instrument handles.

- Pushwires: Manual instrument component of the Biplan instrument set.

- Sterilization trays: Trays designed to hold osteosynthesis screw sets or osteosynthesis screw instrument sets and can serve as sterilization containers.

2. PRELIMINARY WARNINGS

Carefully read and follow the warnings given below. Failure to follow the warnings described in these Ifus will result in the immediate voiding of any warranty on the medical device. BRM Extremities S.r.l. shall not be liable for any damage to persons or property because of failure to follow the instructions for use.

BEFORE USE, VISUALLY CHECK THE INTEGRITY OF THE DEVICES AND IN PARTICULAR THAT THEY SHOW NO SIGNS OF DAMAGE, DEFORMATION OR CORROSION AND THAT THE MARKINGS ARE CORRECTLY LEGIBLE. IF NOT, DO NOT USE AND CONTACT THE MANUFACTURER FOR REPLACEMENT OF THE DEVICES.

7. USE OF THE DEVICE

After being sterilized, the devices can be used for the application of interest. Their handling is restricted to qualified personnel.

DO NOT USE OR REPROCESS OR IN CASE OF POSSIBLE LOSS OF STERILITY, THE DEVICE MUST BE REPROCESSED. IMPLANTATION OF POTENTIALLY CONTAMINATED DEVICES COULD LEAD TO INFECTION, RESULTING IN THE NEED FOR MEDICAL TREATMENT AND POSSIBLE IMPLANT FAILURE.

7.2 METHOD OF USE

The devices are intended to prepare the anatomical area and implant the instruments, as shown in Table 1, according to the surgical technique defined for each device. The indications for use, the contraindications and the instructions for use are summarized below.

HANDLES—Mount the required instrument on the handle, ensure that the coupling is secure and proceed to use the instrument according to surgical technique. The liberation of the handle is performed by removing the instrument from the handle.

THE BLADE HOLDER HANDLES ARE COMPATIBLE WITH BLADE UP TO A MAXIMUM WIDTH OF 3.40 MM AND THICKNESS 0.60 MM.

K-WIRE TUBES—Insert Kirschner wires into the tube prior to sterilization. Once sterilized, as required by the surgical technique, open the tube and remove the required threads.

SLIDING UNIT—Place the bearing unit on the BRDNPHTANDLE handle, beat by sliding the instrument onto the bearing unit, then tap on the end, applying light percussion, remove the sliding unit from the handle.

ADAPTERS—Adapters are used to connect rasps to the available active medical device. They are connected to a rasp via the quick-connect connection and to the active device via

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TRAYS—Insert the intended instruments into the trays, following the indications therein. Close the tray and sterilize it according to the instructions in the Recommended processing instructions for medical devices and accessories intended to be conditioned by the user. Once sterilization is complete, the instruments are ready for use.

3.3 INTENDED PATIENT GROUP

These devices are accessories to the implantable devices shown in Table 1. For identification of the intended patient group, please refer to the Ifu of the reference implantable device.

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These devices are accessories to the implantable devices shown in Table 1. For identification of the indication and expected clinical benefits, please refer to the Ifu of the reference implantable device.

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These devices are accessories to the implantable devices shown in Table 1. For identification of the contraindications, please refer to the Ifu of the reference implantable device.

5.8 CONTOREDIACON

These devices are accessories to the implantable devices shown in Table 1. For identification of the indications and expected clinical benefits, please see the Ifu of the reference implantable device.

5.9 TRATTAMENTO AFTER USE

After use, the devices must be subjected to appropriate reconditioning and control activities, as described in the Recommended processing instructions for medical devices and accessories intended to be conditioned by the user.

10. SHELF LIFE

The devices are supplied in a non-sterile state. The trays are wrapped in a transparent plastic bag, while the remaining devices are packed in transparent polyethylene bags.

TAKE THE DEVICE OUT OF THE PACKAGING BEFORE PROCEEDING WITH REPROCESSING OPERATION.

IF THE PACKAGING IS DAMAGED, DO NOT USE THE DEVICE AND CONTACT THE MANUFACTURER.

6. DEVICE CONDITIONING

Devices should be cleaned and sterilized before use, following instructions in the documentation for connection to an active device and a hollow tip, usually shaped, on which a titanium alloy cylinder slides, which forms a connection system for rasps and quick-connect inserters.

- S-type: type characterised by a cylindrical shank with retention recesses for connection to the surgical motor with reciprocal saw system type "Striker" Rem® tm, CORE; TPS tm and Command 2; the quick-connect connection is anodised in yellow and has a laser marking for rapid identification.

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BRM Extremities S.r.l., info@brm-extremities.com Tel. +39 03411693087