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INSTRUCTIONS FOR USE

1 IDENTIFICATION OF THE DEVICES

These Instructions for Use (IFU), provided in hard copy, accompany Class IIa surgical instruments manufactured by BRM Extremities S.r.l. and intended for implantation of BRM Extremities Siltkoe devices. These devices, listed in Table 1, are considered accessories to the BRM Extremities Siltkoe implantable device and may be used only for activities related to the preparation of the anatomical area and implantation of the indicated device.

For any clarifications or requests, contact the Manufacturer at the following address:

BRM Extremities S.r.l.
Registered Office: Via Lorenzo Mascheroni, 29 – 20145 Milano (MI) - Italy
• Operational Headquarters: Via Privata ai Campi, 7 - 23862 Civate (LC) – Italy
• info@brm-extremities.com
• tel. +39 0341 1693087

1.1 GENERAL INFORMATION

The IFU refer to Class IIa surgical instruments used in the preparation of the anatomical area and implantation of BRM Extremities Siltkoe devices. The devices are supplied in NON-STERILE form and are reusable, after conditioning. The devices can be divided into the following categories:

- rasps;
- trial.

Please refer to the implantable device IFU and surgical technique for a list of all instruments to be used for implanting the device.

2 PRELIMINARY WARNINGS

Carefully read and follow the warnings below. Failure to comply with the warnings described in these IFUs will result in the immediate invalidation of any warranty on the medical device. BRM Extremities Srl is not responsible for any damage to people or things resulting from failure to comply with the warnings.

Before use, read and understand the instructions for use (IFU) entirely. Failure or negligence in complying with the instructions for use could cause the device to malfunction, with possible damage and injury to the user.

Devices must be checked, washed and sterilized before use, as described in section 6.

Follow the instructions for use described in the IFU: uses different than those intended could expose the patient or user to unconsidered dangers.

Do not modify he device in any way.

The devices are intended for use exclusively by orthopaedic surgeons with clinical experience in the implant surgical technique.

Follow the surgical technique indicated by the manufacturer and use instructions only for their intended purpose.

Check the integrity of the device. In case of damage, do not use and contact the manufacturer for replacement. It is strictly forbidden to use these instruments in combination with other devices/accessories not mentioned in these instructions (section 8), as geometric-dimensional compatibility may not be ensured.

3 INTENDED USE

The intended use of the devices is the preparation of the anatomical area and implantation of BRM Extremities Siltkoe device, according to the surgical technique defined by the manufacturer and in combination with the other instrument components.

Refer to the surgical technique for the implantation of the Siltkoe device in order to identify the step of use of these devices.

3.1 INTENDED PATIENT GROUP

These devices are accessories of the Siltkoe implantable device. For identification of the intended patient group, please refer to the IFU of this implantable device.

3.2 INDICATION AND EXPECTED CLINICAL BENEFITS

The surgical instruments in question are accessories to the implantation of the Siltkoe device. It is not possible to define a clinical benefit for an instrument taken individually. These instruments are part, together with others belonging to classes I and/or II, of a set whose clinical benefit consists in the simplification of the related implantation operations. The expected clinical benefit can therefore only be assessed in connection with the implantation of the medical device.

3.3 CONTRAINDICATION

These devices are accessories of the Siltkoe implantable device. For identification of contraindications, please refer to the IFU of the implantable device.

3.4 RESIDUAL RISKS AND SIDE EFFECTS

Given the intended use and the materials used to manufacture the devices, the development of allergic reactions related to the chemical elements present in metal alloys is extremely rare. In the case of particularly sensitive individuals, please contact the manufacturer. No other side-effects related to these devices are known.

In the event of complete loss of function of an instrument during surgery, e.g., due to a fall, it is the surgeon’s responsibility to evaluate the possibility of proceeding with alternative methods or discontinuing the procedure.

The following residual risks related to the individual instruments are also identified:

- Failure to apply the surgical technique prescribed by the manufacturer could result in injury to the patient, including permaent injury.
- In the event of complete loss of function of an instrument during surgery, e.g., due to a fall, it is the surgeon’s responsibility to evaluate the possibility of proceeding with alternative methods or discontinuing the procedure.
- The indicated processing method may not be suitable for complete prion inactivation. Follow local or national guidelines.

4 INFORMATION ON IMPLANTABLE DEVICES

N.A. as the devices in question are not implantable.

5 DEVICE DESCRIPTION AND EXPECTED PERFORMANCE

5.1 RASPS

Rasps are motor-driven steel instruments consisting of a cylindrical shank and a truncated pyramid-shaped end, geometrically and dimensionally similar to that of the relevant implantable device, on which abrasive cusps are placed. These are connected to the handle or,

either directly or via the appropriate adapter, to a reciprocating surgical motor, and are used to remove the bone and obtain a suitable cavity to house the instrument to be implanted.

They are available in two shank versions

- cylindrical shank with retention niches for anchoring to the Multifunctional Siltkoe Handle and for direct anchoring to Striker® Riemtù tm, CORE, TPS tm and Command 2 reciprocating saw systems;
- cylindrical shank with retaining pin for anchoring Multifunctional Siltkoe Quick Handle or to Quick Connection Adapter devices.

For each shank type, they are available in three sizes (20-30, 35, 40-50) for both distal and proximal variants.

5.2 TRIALS

Trial spacers, also called Trial spacers, are coloured medical silicone devices with the same geometry as Siltkoe implantable devices, used by the surgeon to confirm that the chosen device size is correct. They are coloured in different colours according to size and available in 5 right and 5 left sizes.

The following table shows the colour legend.

Trial	Code	Colour
Trial size 20	BRDNPTSTZ0L / BRDNPTSTZ0R	Colour GREEN
Trial size 30	BRDNPTST30L / BRDNPTST30R	Colour ORANGE
Trial size 35	BRDNPTST35L / BRDNPTST35R	Colour YELLOW
Trial size 40	BRDNPTST40L / BRDNPTST40R	Colour BLUE
Trial size 50	BRDNPTST50L / BRDNPTST50R	Colour RED

The colour code can also be found on the label of the Siltkoe implantable device.

5.3 MATERIALS CONSTITUTING THE DEVICE

The rasps are made of AISI 630 stainless steel and the test spacers are made of medical grade silicone.

5.4 PACKAGING AND STERILITY

The devices are supplied non-sterile and individually packed in transparent polyethylene bags.

TAKE THE DEVICE OUT OF THE PACKAGING BEFORE PROCEEDING WITH REPROCESSING OPERATIONS.

IF THE PACKAGING IS DAMAGED, DO NOT USE THE DEVICE AND CONTACT THE MANUFACTURER.

6 DEVICE CONDITIONING

Devices should be cleaned and sterilized before use, following instructions in the document **Recommended processing instructions for medical devices and accessories intended to be conditioned by the user** attached to this IFU. If the conditioning instructions are lost, contact the manufacturer.

7 USE OF THE DEVICE

7.1 PREPARING THE DEVICE

After being sterilized, the device can be used for the application of interest. Their handling is restricted to qualified personnel.

EXERCISE EXTREME CAUTION TO AVOID CONTAMINATION OF THE DEVICE PRIOR TO USE. IF THE DEVICE IS DROPPED OR A LOSS OF STERILITY IS POSSIBLE, THE DEVICE MUST BE REPROCESSED. THE USE OF POTENTIALLY CONTAMINATED DEVICES COULD CAUSE INFECTION, RESULTING IN THE NEED FOR MEDICAL TREATMENT AND POSSIBLE IMPLANT FAILURE.

7.2 MODES OF USE

The devices are intended for the implantation of the Siltkoe device, according to the relevant surgical technique. The operating methods of the individual instruments are summarized below.

BEFORE USE, CHECK THAT THE DEVICES SHOW NO SIGNS OF DAMAGE OR CORROSION AND THAT THE MARKINGS ARE CORRECTLY LEGIBLE. IF THIS IS NOT THE CASE, CONTACT THE MANUFACTURER TO REPLACE THE DEVICES.

RASPS HAVE SHARP ELEMENTS. HANDLE WITH CARE.

IN CASE OF DEVICE BREAKAGE WITH RELEASE OF SMALL FRAGMENTS, REMOVE VISIBLE FRAGMENTS AND VERIFY BY RADIOGRAPHIC TECHNOLOGY THE ABSENCE OF RESIDUAL FRAGMENTS.

7.1 RASPS

- Take the rasp size consistent with the cutting guide used for resection.
- Connect the rasp to the handle or surgical motor, using the adaptors if necessary, and check that the coupling is properly tightened.
- Insert the rasp into the metatarsal and phalangeal medullary canal, taking great care to keep the rasp on axis (keep the upper and lateral and medial edges parallel to the edges of the relevant bone portions).
- Insert the rasp until the stop on the rasp is reached.

BEFORE USING THE RASP, CHECK THAT THE HANDLE, ADAPTOR OR MOTOR IS PROPERLY TIGHTENED.

7.2 TRIALS

- Take the trial size consistent with the size of device you intend to implant. If in doubt between different sizes, always test the smallest size first.
- Insert the trial spacer in the anatomical area and verify the perfect adherence of the trial to the resected bone surfaces, joint mobility and implant stability.

ALTHOUGH THEY ARE GEOMETRICALLY EQUIVALENT TO IMPLANTABLE SILKTOE SPACERS, THE TRIALS ARE MADE OF MATERIAL NOT SUITABLE FOR LONG-TERM IMPLANTATION. IMPLANTATION OF SILICONE TRIAL SPACERS IS STRICTLY PROHIBITED.

8 DEVICES INTENDED TO BE USED TOGETHER

The devices to be used together are summarized in the following table.

Device	Devices intended to be used together		
BRDNPTR2030P BRDNPTR2030D BRDNPTR3535P BRDNPTR3535D BRDNPTR4050P BRDNPTR4050D	20-30 Siltkoe Metatarsal Rasp 20-30 Siltkoe Phalangeal Rasp 35 Siltkoe Metatarsal Rasp 35 Siltkoe Phalangeal Rasp 40-50 Siltkoe Metatarsal Rasp 40-50 Siltkoe Phalangeal Rasp	BRDNPHANDLE Surgical motor	Multifunctional Siltkoe Handle
BRDNPTR2030P-QK BRDNPTR3535P-QK BRDNPTR4050P-QK BRDNPTR2			

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