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INSTRUCTION FOR USE

1 IDENTIFICATION OF THE DEVICES

These Instructions for Use (IFU), provided in hard copy, accompany Class Ia surgical instruments manufactured by BRM Extremities S.r.l. and intended for implantation of BRM Extremities Osteosynthesis Screws. These devices, listed in Table 1, are considered accessories to the BRM Extremities Osteosynthesis Screws implantable device and may be used only for activities related to the preparation of the anatomical area and implantation of the indicated device.

For any clarifications or requests, contact the Manufacturer at the following addresses:

BRM Extremities S.r.l.
Registered Office: Via Lorenzo Mascheroni, 29 - 20145 Milano (MI) - Italy
• Operational Headquarters: Via Privata ai Campi, 7 - 23862 Civate (LC) - Italy
• info@brm-extremities.com
• Tel. +39 0341 1693087

1.1 GENERAL INFORMATION

These IFU refer to Class Ia surgical instruments intended for use in the preparation of the anatomical area and implantation of BRM Extremities Osteosynthesis Screws devices. The devices are supplied NON-STERILE and must be sterilized before use. The tips and gauges are reusable after reprocessing, while the Kirschner wires are disposable.

The devices can be divided into the following categories:

- drill bits;
- depth gauges;
- Kirschner wires.

Please refer to the implantable device IFU and surgical technique for a list of all instruments to be used for preparing the anatomical area and implanting the devices.

2 PRELIMINARY WARNINGS

Carefully read and follow the warnings below. Failure to comply with the warnings described in these IFUs will result in the immediate invalidation of any warranty on the medical device. BRM Extremities Srl is not responsible for any damage to people or things resulting from failure to comply with the warnings.

Before use, read and understand the instructions for use (IFU) entirely. Failure or negligence in complying with the instructions for use could cause the device to malfunction, with possible damage and injury to the user.

Devices must be checked, washed and sterilized before use, as described in section 6.

Follow the instructions for use described in the IFU: uses different than those intended could expose the patient or user to unconsidered dangers.

Do not modify the device in any way.

The devices are intended for use exclusively by orthopaedic surgeons with clinical experience in the implant surgical technique.

Follow the surgical technique indicated by the manufacturer and use instruments only for their intended purpose.

Check the integrity of the device. In case of damage, do not use and contact the manufacturer for replacement.

It is strictly forbidden to use these instruments in combination with other devices/accessories not mentioned in these instructions (section 8), as geometric-dimensional compatibility may not be ensured.

3 INTENDED USE

The intended use of the devices is the preparation of the anatomical area and implantation of BRM Extremities osteosynthesis screws, according to the surgical technique defined by the manufacturer and in combination with the other instrument components.

Refer to the surgical technique for the implantation of the device Osteosynthesis screws in order to identify the step of use of these devices.

3.1 INTENDED PATIENT GROUP

These devices are accessories of the implantable device Osteosynthesis Screws. For identification of the intended patient group, please refer to the IFU of this implantable device.

3.2 INDICATION AND EXPECTED CLINICAL BENEFITS

The surgical instruments in question are accessories to the implantation of the osteosynthesis screw device. It is not possible to define a clinical benefit for an instrument taken individually. These instruments are part, together with others belonging to classes I and/or I, of a set whose clinical benefit consists in the simplification of the related implantation and/or explanation operations. The expected clinical benefit can therefore only be assessed in connection with the implantation of the medical device.

3.3 CONTRAINDICATION

These devices are accessories of the implantable device Osteosynthesis Screws. For identification of contraindications, please refer to the IFU of the implantable device.

3.4 RESIDUAL RISKS AND SIDE EFFECTS

Given the intended use and the materials used to manufacture the devices, the development of allergic reactions related to the chemical elements present in metal alloys is extremely rare. In the case of particularly sensitive individuals, please contact the manufacturer. No other side effects related to these devices are known.

Intra-operative complications associated with the surgical procedure as a whole, of which the individual instruments in this IFU are an integral part, are reported in the IFU of the individual implantable device.

The following residual risks related to the individual instruments are also identified:

- Failure to apply the surgical technique prescribed by the manufacturer could result in injury to the patient, including permanent injury.
- In the event of complete loss of function of an instrument during surgery, e.g., due to a fall, it is the surgeon's responsibility to evaluate the possibility of proceeding with alternative methods or discontinuing the procedure.
- The indicated reprocessing method may not be suitable for complete prion inactivation. Follow local or national guidelines.

4 INFORMATION ON IMPLANTABLE DEVICES

N.A. as the devices in question are not implantable.

5 DEVICE DESCRIPTION AND EXPECTED PERFORMANCE

5.1 DRILL BITS

The drills are cannulated cylindrical steel instruments; they are equipped with a quick-connect coupling at one end for connection to a handle or surgical motor, and with a tip with cutting elements at the other. They can be used either manually

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Tabelle 1 / Table 1 / Tableau 1 / Tabela 1/ Taha 1/ Tabelle 1

Typologia / Type of Instrument / Tipologia / Tipologia / Typologie	Codice / Code / Código / Código / Code	Denominazione commerciale / Commercial name / Nom commercial / Nome comercial / Handelsname
	STOI217000P	Canulated Drill Bit Ø 1.7mm KW 10
	STOI217008P	Canulated Drill Bit Ø 1.7mm KW 08
	STOI222000P	Canulated Drill Bit Ø 2.2mm
	STOI227012P	Canulated Drill Bit Ø 2.7mm KW 12
	STOI227000P	Canulated Drill Bit Ø 2.7mm KW 15
	STAI239000P	Canulated Drill Bit Ø 3.9mm
	STAI246000P	Canulated Drill Bit Ø 4.6mm
	STAI255000P	Canulated Drill Bit Ø 5.5mm
	STVI220000P	Canulated Drill Bit Ø 1.8mm
	STVI224000P	Canulated Drill Bit Ø 2.0mm
	STVI234000P	Canulated Drill Bit Ø 2.9mm
	STVI220XLOP	Long Canulated Drill Bit Ø 1.8mm
	STVI224XLOP	Long Canulated Drill Bit Ø 2.0mm
	STVI234XLOP	Long Canulated Drill Bit Ø 2.9mm
	STAI239GDDP	Canulated Drill Bit Ø 3.9mm for Guidedriß
	STAI246GDDP	Canulated Drill Bit Ø 4.6mm for Guidedriß
	STOI222000M	Depth Gauge 40 mm
	STOI223000M	Depth Gauge 60 mm
	STAI256000M	Depth Gauge 85 mm KW 120 mm
	STAI256150M	Depth Gauge 85 mm KW 150 mm
	STAI270000M	Depth Gauge 110 mm
	STOKW08800003N	K-Wire Ø 0.8 x 80 mm
	STOKW09800003N	K-Wire Ø 0.9 x 80 mm
	STOKW10800003N	K-Wire Ø 1.0 x 80 mm
	STOKW08120003N	K-Wire Ø 0.8 x 120 mm
	STOKW08120803N	K-Wire Ø 0.8 x 120 (80mm Marker)
	STOKW09120003N	K-Wire Ø 0.9 x 120 mm
	STOKW09120803N	K-Wire Ø 0.9 x 120 (80mm Marker)
	STOKW10120003N	K-Wire Ø 1.0 x 120 (80mm Marker)
	STOKW10150003N	K-Wire Ø 1.0 x 150 (80mm Marker)
	STOKW12080003N	K-Wire Ø 1.2 x 80 mm
	STOKW12120003N	K-Wire Ø 1.2 x 120 mm
	STOKW12120803N	K-Wire Ø 1.2 x 120 (80mm Marker)
	STOKW12150003N	K-Wire Ø 1.2 x 150 mm
	STOKW12150803N	K-Wire Ø 1.2 x 150 (80mm Marker)
	STOKW15080003N	K-Wire Ø 1.5 x 80 mm
	STOKW15120003N	K-Wire Ø 1.5 x 120 mm
	STOKW15150003N	K-Wire Ø 1.5 x 150 mm
	STOKW16150803N	K-Wire Ø 1.6 x 150 (80mm Marker)
	STOKW16120003N	K-Wire Ø 1.6 x 120
	STOKW16150003N	K-Wire Ø 1.6 x 150 mm
	STOKW16151203N	K-Wire Ø 1.6 x 150 (120mm Marker)
	STOKW16200003N	K-Wire Ø 1.6 x 200 (120mm Marker)
	STOKW16201503N	K-Wire Ø 1.6 x 200 (150mm Marker)
	STOKW20200003N	K-Wire Ø 2.0 x 200 mm
	STOKW25200003N	K-Wire Ø 2.5 x 200 mm
	STOKW08800010N	K-Wire Ø 0.8 x 80 mm
	STOKW09800010N	K-Wire Ø 0.9 x 80 mm
	STOKW10800010N	K-Wire Ø 1.0 x 80 mm

	STOKW08120010N	K-Wire Ø 0.8 x 120 mm
	STOKW08120810N	K-Wire Ø 0.8 x 120 (80mm Marker)
	STOKW09120010N	K-Wire Ø 0.9 x 120 mm
	STOKW09120810N	K-Wire Ø 0.9 x 120 (80mm Marker)
	STOKW10120010N	K-Wire Ø 1.0 x 120 mm
	STOKW10120810N	K-Wire Ø 1.0 x 120 (80mm Marker)
	STOKW10150010N	K-Wire Ø 1.0 x 150 mm
	STOKW10150810N	K-Wire Ø 1.0 x 150 (80mm Marker)
	STOKW12080010N	K-Wire Ø 1.2 x 80 mm
	STOKW12120010N	K-Wire Ø 1.2 x 120 mm
	STOKW12150010N	K-Wire Ø 1.2 x 150 (80mm Marker)
	STOKW12150810N	K-Wire Ø 1.2 x 150 (80mm Marker)
	STOKW15080010N	K-Wire Ø 1.5 x 80 mm
	STOKW15120010N	K-Wire Ø 1.5 x 120 mm
	STOKW15120810N	K-Wire Ø 1.5 x 120 (80mm Marker)
	STOKW15150010N	K-Wire Ø 1.5 x 150 mm
	STOKW15150810N	K-Wire Ø 1.5 x 150 (80mm Marker)
	STOKW16120010N	K-Wire Ø 1.6 x 120
	STOKW16150010N	K-Wire Ø 1.6 x 150 mm
	STOKW16151210N	K-Wire Ø 1.6 x 150 (120mm Marker)
	STOKW16200010N	K-Wire Ø 1.6 x 200 mm
	STOKW16201510N	K-Wire Ø 1.6 x 200 (150mm Marker)
	STOKW20200010N	K-Wire Ø 2.0 x 200 mm
	STOKW25200010N	K-Wire Ø 2.5 x 200