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		Countersink for FX screws with Ø 2,2-2,6
		BTX7 Screwdriver (solid or cannulated)
		Wire guide 1,2 mm
		Parallel Guidewire
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

		Tip Ø 2,2 mm
		Countersink for screws with Ø 3,2
		BTX8 Screwdriver (solid or cannulated)
		Wire guide 1,2 mm
		Parallel Guidewire
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

SOS132BBBYW	Cannulated osteosynthesis screws TOOL - Ø 3,2	Tips Ø 2,7 mm
		Countersink for screws with Ø 4
		BTX9 Screwdriver (solid or cannulated)
		Wire guide 1,2 mm
		Parallel Guidewire
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

SOS140BBBYW	Cannulated osteosynthesis screws TOOL - Ø 4	Tips Ø 3,9 mm
		Countersink for screws with Ø 5,6
		BTX20 Screwdriver (solid or cannulated)
		Wire guide 1,6 mm
		Depth Gauge for reference at 120 mm and 150
		Kirschner Wires (according to Table 5)

SOS156BBBYW	Cannulated osteosynthesis screws TOOL PRO Ø 5,6	Tips Ø 4,6 mm and Ø 5,5 mm
		Countersink for screws with Ø 7
		BTX25 Screwdriver (solid or cannulated)
		Wire guide 2,5 mm
		Tip guide 4,6 mm
		Tissue protection sleeve TOOL PRO Ø 7,0
		Depth Gauge for reference at 200 mm
		Kirschner Wires (according to Table 5)

SOS170BBBYW	Cannulated osteosynthesis screws TOOL PRO Ø 7	Tips Ø 1,8 mm
		Countersink for TOOL VIP Ø 2,5
		Screwdriver Hex. 1,7 (solid or cannulated)
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

SSL125BBBW	Cannulated osteosynthesis screws TOOL VIP Ø 2,5	Tips Ø 2,9mm
		Countersink for TOOL VIP Ø 4
		Screwdriver Hex. 2,5 (solid or cannulated)
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

SSL130BBBW	Cannulated osteosynthesis screws TOOL VIP Ø 3	Tips Ø 2mm
		Countersink for TOOL VIP Ø 3
		Screwdriver Hex. 2 (solid or cannulated)
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

SSL140BBBW	Cannulated osteosynthesis screws TOOL VIP Ø 4	Tips Ø 2,9mm
		Countersink for TOOL VIP Ø 4
		Screwdriver Hex. 2,5 (solid or cannulated)
		Depth Gauge for reference at 80 mm
		Kirschner Wires (according to Table 5)

SROXAABBBW	VR Screws	VR Screws
		Screwdriver for VR Screws
		In addition, there are three instruments, not specific to one device rather than another, which are useful for preparing the implant area: double rasp, straight and curved elevator.
		The surgical instruments are preferentially supplied in the form of a set, so that the equipment required for the relevant surgical technique is immediately available; however, they can also be supplied individually. In this case, it is recommended to contact the manufacturer for a list of the individual instruments to be used together or for clarification.
		Specific instrument sets have been developed per screw type and diameter. They are summarised in Table 4.

Table 4 – Devices intended to be used together	DEVICE	DEVICES INTENDED TO BE USED TOGETHER
	STO122BBYS	Cannulated osteosynthesis screws TOOL - Ø 2,2-2,6
	SET.TOOL2240P	Instrument box for screws Ø 2,2-4,0 small
	STO132BBYS	Cannulated osteosynthesis screws TOOL - Ø 3,2
	SET.TOOL2240P	Instrument box for screws Ø 2,2-4,0 small
	STO140BBYS	Cannulated osteosynthesis screws TOOL - Ø 4
	SET.TOOL2240P	Instrument box for screws Ø 2,2-4,0 small
	SOS156BBBYW	Cannulated osteosynthesis screws TOOL PRO
	SET.TOOL5670P	Instrument box for screws TOOL PRO
	SSL125BBBW	Cannulated osteosynthesis screws TOOL VIP Ø 2,5
	SET.TOOLVIPP	Instrument box for screws TOOL VIP
	SSL130BBBW	Cannulated osteosynthesis screws TOOL VIP Ø 3
	SET.TOOLVIP	Instrument box for screws TOOL VIP with screw holder
	SFX122BBBS	Cannulated osteosynthesis screws TOOL FX
	SET.TOOLXP	Instrument box for screws TOOL FX
	SRD120BBBS	VR Screws
	SET.VR	Instrument box for break-off screws

In addition, there are three instruments, not specific to one device rather than another, which are useful for preparing the implant area: double rasp, straight and curved elevator.

The surgical instruments are preferentially supplied in the form of a set, so that the equipment required for the relevant surgical technique is immediately available; however, they can also be supplied individually. In this case, it is recommended to contact the manufacturer for a list of the individual instruments to be used together or for clarification.

Specific instrument sets have been developed per screw type and diameter. They are summarised in Table 4.

For the correct execution of the operation, the surgical technique also involves the use of Kirschner wires combined according to the scheme shown in Table 5

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ATTENTION: USE KIRSCHNER WIRES COMPATIBLE WITH THE SCREW TO BE IMPLANTED, ACCORDING TO TABLE 5

9 TREATMENT AFTER USE

The device is for single use only and is supplied in sterile packaging. If the packaging is opened and the sterile barrier is lost, the device must be disposed of.

10 SHELF LIFE

The shelf life of the device is 5 years from the date of sterilisation. The expiry date and date of manufacture are indicated on the pouch label.

11 DISPOSAL OF DEVICE AND PACKAGING

When disposing of the device, follow the waste management regulations of the local authority (special hospital waste). Do not under any circumstances dispose of the product or its packaging in the environment.

Packaging

Materal family

Decision 97/129/CE

Envelope Plastic - Aluminium C/LDPE90 >OPA+PE-7

Box Non-corrigated cardboard PAP21

IFU Paper PAP22

Protective film Polyethylene LDPE4

12 TRANSPORT AND STORAGE CONDITIONS

The devices can be transported and stored at room temperature and must be kept dry, away from heat sources and away from light.

13 INCIDENT REPORTING AND SURVEILLANCE ACTIVITIES

WARNING: In accordance with current legislation, the user is obliged to report any serious incidents occurring with the device to the

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