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RECOMMENDED PROCESSING INSTRUCTIONS FOR MEDICAL DEVICES AND ACCESSORIES

INTENDED TO BE CONDITIONED BY THE USER Instructions in accordance with ISO 17664 and AAMI S781. **Manufacturer:** BRM Extremities S.r.l. **Method:** Steam sterilization (moist heat)

Application field: These instructions apply to:

- single use medical devices supplied by BRM Extremities in a non-sterile version and intended to be sterilized by the user before use;

- reusable surgical instruments supplied by BRM Extremities and intended for reprocessing in a reusable facility.
- Unless explicitly noted in the instructions accompanying a specific instrument, all BRM Extremities devices and accessories may be safely and effectively processed according to the cleaning instructions (manual or manual / automated mode) and the sterilization parameters provided in this document.

In countries where the reprocessing requirements are more stringent than those provided in this document, the user / processor will be responsible for complying with the most stringent legislation.

BRM Extremities validated the reprocessing of surgical devices and instruments but, in order to obtain the desired result, it is the responsibility of the user / hospital / health care worker to ensure that the conditioning takes place using suitable equipment and materials and that the staff has been adequately trained; this requires equipment and processes to be validated and checked regularly. Any deviation from these instructions by the user / hospital / healthcare staff will have to be evaluated for effectiveness to avoid potential adverse consequences.

WARNINGS
BRM Extremities devices provided in a NON-STERILE version must be cleaned and sterilized according to these instruction prior to use.

- The devices Osteosynthesis Screws and Kirschner Wires are single use devices. Once they have come into contact with the patient, they cannot be reprocessed but must be disposed of according to current regulations.

- Any protective packaging must be removed from the devices before the first cleaning and sterilization cycle.

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Ethylene Oxide (EO), gas plasma and dry heat sterilization methods are not recommended for sterilization of BRM Extremities medical devices. Steam (moist heat) is the recommended method.

- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated devices.

- Caution should be exercised while handling, cleaning or wiping instruments with cutting edges, tips and teeth.

- Safety cleaning disinfection agents containing aldehyde, chloride, active chlorine, bromine, bromide, iodine or iodine are corrosive and should not be used.

- Prevent any biological impurities from drying on contaminated devices. All subsequent cleaning and sterilization steps are facilitated by preventing blood, body fluids and tissue debris from drying on used instruments.

- Metal brushes and scouring pads must not be used during manual cleaning. These materials will damage the surface and finish of the devices. Use only soft bristle nylon brushes with different shapes, lengths and sizes to aid with manual cleaning.

- When processing instruments with polymer components, avoid the use of plastic instruments.

- Use of hard water should be avoided. Softened tap water may be used for the preparation of detergent solutions and intermediate rinsing, however critical water should be used for final rinsing to prevent mineral deposits.

- For reprocessing activities, only use cleaning agents that are classified as medical devices in accordance with EU Regulation 2017/745 or Directive 93/42/CE and have been certified accordingly. The use of enzymatic proteolytic detergents is recommended.

Do not process instruments with polymer components at temperatures equal to or greater than 140 °C/285°F because severe surface damage to the polymer may occur.

- Oils or silicone lubricants should not be used on medical devices.

- All surgical instruments made of steel contain nickel in alloy form: users with proven evidence of contact dermatitis due to allergic reaction to nickel are recommended to use personal protective equipment (gloves).

I) PACKAGING FOR STERILIZATION

Single devices may be packaged in a medical grade sterilization pouch or wrap compliant with ISO 11607. **Please note that only FDA cleared pouches/wraps are recommended for use in the US.** During packing, check the pouch or sheet for tears. The devices must be double-wrapped or by an equivalent method (ref.: AAMI S779, AORN guidelines). The use of reusable packaging is not recommended.

The set of devices must be placed in the dedicated tray (accessory of the medical device) by inserting each individual instrument in the specific appropriately identified housing. The box or tray must be wrapped in medical grade, FDA approved and / or ISO 11607 compliant sterilization wrap following the double wrap method or equivalent method (ref.: AAMI S779, AORN guidelines). **Please note that only FDA cleared pouches/wraps are recommended for use in the US.**

L) STERILIZATION
Moist/steam sterilization is the recommended method for sterilization. The use of an approved chemical indicator (type 5) or chemical emulato (type 6) within each sterilization load is recommended. For loading configuration and operation of the sterilisation equipment and its maintenance, always consult and follow the manufacturer's instructions.

NOTE: DO NOT STACK TRAYS DURING STERILIZATION

Notwithstanding the above, the devices may be sterilized in autoclave (e.g. FDA clearance, EN 13060 or EN 285 compliance) and must be classified as a medical device in accordance with EU Regulation 2017/745 or Directive 93/42/CE. **Please note that only FDA cleared devices are recommended for use in US markets.** The effectiveness of sterilization is validated according to the following exposure times and temperatures, aimed at guaranteeing the achievement of a SAL sterility level of 10⁻⁶.

Gravity Cycle: 121°C—exposure time 30 minutes—drying time: 15-30 min.

Vacuum pulse (not for U.S.A. market): 134°C—exposure time 4 min—drying time: 20 min.

Both methods have been validated by the manufacturer; the choice of cycle must be made by the user, based on good practice and available equipment. The sterilisation process adopted must always be validated by the user, in accordance with the EN ISO 17665-1. The use of sterilization parameters or methods other than those indicated is the sole responsibility of the user, who must implement appropriate checks to ensure the effectiveness of the sterilization and the maintenance of the mechanical and dimensional characteristics of the sterilized product.

Drying and cooling

The recommended drying time for single wrapped instruments is 20 minutes unless otherwise stated in device specific instructions.

Drying times for instruments processed in wrapped trays can vary depending upon the type of packing, type of instruments, type of sterilizer and total load. A minimum dry time of 30 minutes is recommended, but to avoid wet packs, extended dry times greater than 30 minutes may be needed for larger loads under certain conditions or if otherwise recommended in accompanying documentation. For large loads verification of dry times by the health care provider is recommended.

A 30 minutes minimum cooling time is recommended after drying but longer times may be used. Warning: Instruments should be cooled in a basin containing critical (reverse osmosis or distilled) water for at least two minutes or until all traces of detergent or biological residues have disappeared. During rinsing, operate all moving and hinged parts. Rinse thoroughly and vigorously with the help of a syringe, the lumens, holes, cannulas and other areas that are difficult to access.

Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

Step 1: remove all devices from the enzyme solution and rinse with running water for at least two minutes. During rinsing, operate all moving and hinged parts. Thoroughly and vigorously rinse lumens, holes, cannulas and other hard-to-reach areas.

Step 2: completely submerge instruments in the enzyme solution and repeatedly actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulas should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. Keep the devices immersed in the solution for the time indicated by the manufacturer of the detergent and, in any case, not less than 5 minutes.

Step 3: during soaking time, brush the surfaces using a soft nylon bristle brush until all visible dirt is removed. Activate the mobile mechanisms. Pay particular attention to lumen, hinged parts, locks, teeth, surfaces and rough areas with moving parts. Lumens, blind holes and cannulas must be washed with a cylindrical brush with nylon bristles of the appropriate size and with syringes of adequate size, taking the detergent solution directly from the immersion bath.

Note: Do not immerse them in a basin containing critical (reverse osmosis or distilled) water for at least two minutes or until all traces of detergent or biological residues have disappeared. During rinsing, operate all moving and hinged parts. Rinse thoroughly and vigorously with the help of a syringe, the lumens, holes, cannulas and other areas that are difficult to access.

Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

Step 1: remove all devices from the enzyme solution and rinse with running water for at least two minutes. During rinsing, operate all moving and hinged parts. Thoroughly and vigorously rinse lumens, holes, cannulas and other hard-to-reach areas.

Step 2: completely submerge instruments in the enzyme solution and repeatedly actuate instruments with hinges or moving parts to ensure contact of the solution

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