Cod 100025 rev 9 del 25/03/2024

FRESE PER CHIRURGIA PERCUTANEA DEL PIEDE

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BURRS FOR PERCUTANEOUS FOOT SURGERY
USE AND MAINTENANCE MANUAL
This document is the exclusive property of BRM Extremities S.r.l. All rights reserved. Any disclosure, reproduction or transfer of the content to third parties is prohibited without the Company's prior consent.
BRM Extremities Srl declines any responsibility for damage to people or things due to improper use of this product and failure to comply with the indications, warnings, instructions and precautions contained in this user manual.
This manual of use is provided only in paper format and must always go with the medical device.
The burrs for percutaneous foot surgery are manufactured by:

BRM Extremities S.r.l.
Registered office: Via Lorenzo Mascheroni 29 - 20145 Milano (MI) - Italy
Headquarter: Via Privata ai Campi 7 - 23862 Civate (LC) - Italy
P. IVA - 08683610961 Tel. - 0341.169387

1. GENERAL INFORMATION
BRM Extremities S.r.l. burrs are medical devices intended to be used in mini-invasive percutaneous foot surgery. Burrs are available in different head shape and dimensions, both in sterile/disposable and not sterile/reusable (up to 10 uses) version. Burrs are made up in thermally treated AISI 630 stainless steel which guarantees a proper biocompatibility according to the intended use.

1.1 INTENDED USE
The intended use of the burrs is the cutting and/or the milling of the bone tissue (osteotomy), through an exclusively mechanical action, during mini-invasive percutaneous foot surgery.

1.2 CLASSIFICATION
Classification according to the Directive 93/42/CEE Annex IX rule 6 class Ia.

1.5 SYMBOLOLOGY
Below are shown the symbols used for the management of the main warnings for a proper and safe use of the device, in order to make more comfortable and easier the manual reading.

Requirement for a correct use
This symbol identifies the presence of information for correct use of the device.
Information requirement
This symbol identifies the presence of useful and general information whose reading guides the user to a conscious use of the device and / or to the execution of actions.
This symbol identifies that the device is design and manufactured according to the safe requirements (RES) of the Medical Devices Directive 93/42/CEE (Medical device class Ia, according to classification rule n.6 of the Annex IX).

1.6 PRELIMINARY WARNINGS
Failure to comply with the warnings listed below as well as with the standards and precautions described in this user manual entails the immediate termination of any warranty on the Burr's Device. BRM Extremities S.r.l. is not responsible for any damage to people or things as a result of failure to comply with the rules or precautions listed below and reported in this user manual.

LACKS OR NEGLIGENCE IN COMPLIANCE WITH THE FOLLOWING INDICATIONS MAY CAUSE DEVICE MALFUNCTION, DAMAGE AND INJURY TO THE USER OR PATIENT
DO NOT USE THE DEVICE BEFORE THE READING AND UNDERSTANDING OF THIS USE MANUAL IN ALL ITS PARTS
NO MODIFICATION OF THE DEVICE AND/OR OF ITS PARTS IS ALLOWED
THE USE OF THE DEVICE FOR PURPOSES OTHER THAN THOSE INDICATED IN THIS USER MANUAL COULD EXPOSE THE PATIENT AND THE OPERATOR TO HAZARDS

1.7 PACK CONTENTS

Burrs are available in the following configurations:

- Single pack (1 burr) sterile and disposable
- Multiple pack (6 burrs) sterile and disposable
- Single pack (1 burr) not sterile and reusable
- Multiple pack (6 burrs) not sterile and reusable

Regardless of the end configuration, burrs are always individually packed in peelable PET12 / AL9 / PES0 bags and placed in a cardboard box. The secondary packaging always contains this user manual and a series of adhesive labels to be applied on the patient's medical record or on other documents in order to ensure the traceability of the production batch.

In the event that all the elements listed above are not included in the pack, contact the manufacturer immediately.

2. METHOD OF USE

2.1 BURRS TYPES AND THEIR FUNCTIONS

The following kinds of burrs are available:

- Shannon Standard: burr with helical cutting edges (15°), a diameter of 2 mm and a cutting-edge length of 12 mm.
- Shannon Short: burr with helical cutting edges (8°), a diameter of 2 mm and a cutting-edge length of 8,3 mm.
- Shannon Isham: burr with straight cutting edges, a diameter of 2 mm and a cutting-edge length of 12 mm.
- Conical: burr with helical cutting edges (10°), a diameter of 3,1 or 4,1 mm and a cutting-edge length of 13 mm.
- Straight Cylindrical: burr with straight cutting edges, a diameter of 2, 2,5 or 3,1 mm and a relative cutting-edge length of 14, 20 or 22 mm.
- Wedge: burr with helical cutting edges (10°), a diameter of 3,1 or 4,1 mm and a cutting-edge length of 13 mm.
- Rounded Cylindrical: burr with straight cutting edges, a diameter of 3,1 or 4,1 mm and a cutting-edge length of 15 mm.
- Brophy: burr with helical cutting edges (13°), a diameter of 5 mm and a cutting-edge length of 14,8 mm.

Type	Indication
Short and Standard Shannon	Osteotomy of the minor finger phalanges; Exostectomy of the minor fingers; Remodeling of bone surfaces; Dorsolateral wedge-metatarsal exostectomy (standard length).
Shannon Isham	Osteotomy of the minor finger phalanges and of the first metatarsal (distal e proximal); Akin Osteotomy; Exostectomy of the minor fingers; Remodeling of bone surfaces; Calcaneal spur exostectomy; Dorsolateral exostectomy of the head of the fifth metatarsal; Dorsolateral wedge-metatarsal exostectomy.
Wedge	Husking; Bunionectomy of the first metatarsal; Calcaneal spur exostectomy; Dorsolateral wedge-metatarsal exostectomy; Akin Osteotomy.
Straight Cylindrical	Bunionectomy of the first metatarsal; Calcaneal spur exostectomy; Dorsolateral wedge-metatarsal exostectomy.
Rounded Cylindrical	Moderate and severe exostectomy of the big toe
Conical	Cheilectomy (the shape allows to easily remove the osteophytes of the metatarsophalangeal joint)
Brophy	Exostectomy of large surfaces

THE CONTENTS PROVIDED IN THIS CHAPTER ARE INTENDED ONLY TO BE INDICATIVE. THE CHOICE OF SURGICAL TECHNIQUE AND APPROPRIATE INSTRUMENTS FOR EACH PATIENT IS UNDER THE RESPONSIBILITY OF THE SURGEON, IN ACCORDANCE WITH HIS PRACTICE, EXPERIENCE, TRAINING, LEVEL OF ACCURACY AND KNOWLEDGE OF MEDICAL LITERATURE. THE PATIENT SELECTION CRITERIA AND THE ASSESSMENT OF THE USABILITY OF BRM EXTREMITIES BURRS WITHIN THE SPECIFIC CASE IS THE RESPONSIBILITY OF THE SURGEON, AFTER AN ADEQUATE ANAMNESIS, EVEN THE INFORMATION TO THE PATIENT ON THE POSSIBLE RISKS ARISING FROM THE USE OF THE MEDICAL DEVICES OBJECT OF THIS INSTRUCTIONS FOR USE IS RESPONSIBILITY OF THE SURGEON, WITH REFERENCE TO THE IMPLEMENTATION OF THE REGULATIONS CONCERNING THE INFORMED CONSENT.

2.2 EXTRACTION FROM PACKAGING

CHECK THE INTEGRITY OF THE DEVICE AND OF THE PRIMARY PACKAGING BEFORE PROCEEDING WITH THE OTHER STEPS. IF THE DEVICE OR THE PRIMARY PACKAGING ARE DAMAGED DO NOT USE THE DEVICE AND DISPOSE OF IT ACCORDING TO THE INSTRUCTIONS OF THE "DISPOSAL" CHAPTER.
CHECK THE EXPIRATION DATE OVER THE PACKAGING OF THE STERILE BURRS: IF THE DATE HAS BEEN EXCEEDED DO NOT USE THE DEVICE AND DISPOSE OF IT FOLLOWING THE INSTRUCTIONS OF THE "DISPOSAL" CHAPTER.

The handling of the burrs is reserved to qualified personnel (instrumentalist or qualified professional nurse). It is recommended to remove the burrs from their primary packaging only immediately before their use and as far as possible to handle the burrs wearing talc-free gloves in order to not contaminate them.

2.3 BURRS PREPARATION BEFORE USE

2.3.1 Sterile Burrs
Sterile and disposable burrs, when removed from their primary packaging, can be immediately used for the intended use.
2.3.2 Not Sterile Burrs
Not Sterile and reusable burrs must be washed, disinfected and sterilized before their use, according to the following instructions:

- Immerse the burrs in a basin containing warm water and neutral enzymatic detergent, for a minimum time of 5 minutes. This step can be carried out in an ultrasonic tank.
- Note: It is necessary to use a cleaning product suitable for the disinfection of surgical medical devices and to follow the manufacturer's instructions for concentration and contact times, which in any case must not be less than 5 minutes.

DURING WASHING DO NOT USE HYDROGEN PEROXIDE (OXYGEN WATER) THAT COULD DAMAGE THE DEVICES.
2)Place the burrs in a suitable container (steel or Teflon tray or basket) and rinse with sterile water for injections for at least 2 minutes.
3)Sterilization in Autoclave: the recommended method is to package the items or put them in a suitable surgical tray and sterilize for 5 minutes at 134° C and 2,1 bar.

4)Place the sterilized surgical instruments in a tray containing sterile physiological solution during the surgery.
2.4 BURRS USE
To be used the burrs have to be inserted into the mandrel of a surgical micromotor (handpiece). Burr's stem has been designed in conformity with the EN 1797 technical standard: it is therefore recommended to use a handpiece with a compatible mandrel, making sure that the Shank of the burr is well inserted and locked in the mandrel. During burrs use, it is necessary to keep the rotation speed below 8000 rpm in order to avoid overheating of the bone tissue with possible consequent damage (tissue necrosis). To limit overheating, other precautions should also be taken, such as spraying saline at the intervention site.

2.5 BURRS TREATMENT AFTER USE

2.5.1 Sterile Burrs
Sterile burrs are disposable therefore they have to be disposed after the use.
DO NOT RISTERILIZE AND DO NOT REUSE THE DISPOSABLE STERILE BURRS. IMPROPER REUSE MAY CAUSE THE TRANSMISSION OF SERIOUS INFECTIONS.

2.5.2 Not Sterile Burrs
Non sterile burrs are reusable up to 10 times. After their use follow the steps below:

- Rinse the burrs under warm running water for two minutes.
- Immerse the burrs in a basin containing warm water and neutral enzymatic detergent, for a minimum of 5 minutes. This step can be carried out in an ultrasonic tank.

Note: It is necessary to use a cleaning product suitable for the disinfection of surgical medical devices and to follow the manufacturer's instructions for concentration and contact times, which in any case must not be less than 5 minutes.

DURING WASHING DO NOT USE HYDROGEN PEROXIDE (OXYGEN WATER) THAT COULD DAMAGE THE DEVICES.
3) Clean the devices with a soft bristle brush.
4) Rinse with warm water for two minutes.
5) Disinfect the instruments with a sterile 70% ethyl alcohol solution.
6) Allow the solution to evaporate in a protected environment and place the burrs in the sterilization box.
7) Store the box in a dry and clean place, away from sources of contamination.

BEFORE USE AGAIN THE BURRS REPEAT THE STEPS FROM 1 TO 4 DESCRIBED IN THE CHAPTER 2.3.2

3. USEFUL LIFE

The useful life of the devices depends on the kind of burrs: the shelf life of sterile burrs is defined in 60 months (with reference to primary packaging validation and the maintenance of the sterility). Please refer to the expiration date shown on the package to determine the shelf life of the burrs.

Furthermore, sterile burrs are disposable; therefore, they have to be disposed of after their use. Non sterile burrs can be used up to 10 times, corresponding to the max number of sterilization cycles they can be performed.

4. DISPOSAL

The devices have to be disposed according to the waste management regulation of the Local Management Authority (Special hospital waste). In any case do not disperse the product or its packaging in the environment.

Envelope: Polietilene — C/DPPE81
Box: Non corrugated cardboard — C/PAP21
Protective film: Polyethylene — LDPE4

5. ENVIRONMENTAL CONDITIONS
In order to ensure the maintenance of sterility and integrity of the device, the following environmental conditions should be respected.

Conditions	Humidity:	Storage / transport	5 / 95% Ur without condensation
	Atmospheric pressure:	Use	700 at 1020 hPa
		Storage / transport	500 at 1060 hPa

6. SYMBOLS USED ON THE LABEL

Identification of the Manufacturer	
Product code	
Batch number	
Attention, refer to the user manual. Security information	
Caution, to be sold only on prescription from an authorized health care professional	

STERILE R	Sterilized with ionizing radiation
	Do not resterilize
	Do not use if the packaging is damaged
	Do not reuse
	Consult instructions for use
	Use by the expiration date
	Not sterile
	Keep away from light
	Keep dry
	Environmental storage temperature
	Latex free
CE 0426	Compliant with Council Directive 93/42/EEC and subsequent amendments, DM Class Ia. The 4-digit code indicates the device certification Notified Body.

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