DIGITALIS

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BRM Extremities S.r.l. disclaims any liability for damage to persons or objects, due to improper use of this product and failure to comply with the instructions, failure to observe the instructions, warnings, instructions and precautions contained in this user manual.

This user manual is provided printed only and must always accompany the Digitalis device.

Digitalis is manufactured by:

BRM Extremities S.r.l.
Headquarters - Via Lorenzo Mascheroni 29, 20145 Milan (MI) - Italy
Operational Headquarters - Via Privata ai Campi 7, 23862 Civate (LC) - Italy
VAT Number - 08683610961
Tel. - 0341 1693087

1 GENERAL INFORMATION

Digitalis is a medical device designed and manufactured to offer the surgeon specialized in hand surgery a solution for patients suffering from severe rheumatic diseases (arthritis diseases), for whom arthroplasty surgery of the metacarpophalangeal (MCP) or proximal interphalangeal joint (PIP) of the fingers is required. The device consists of a medical grade silicone elastomer spacer (NuSiTM MED-4735), which replaces the joint and consists of two stems.

the distal stem that guarantees anatomical adaptation in relation to the intramedullary canal of the phalanx;

the proximal stem that helps in joining to the metacarpal (in case of MCP joint) or phalangeal (in case of PIP joint) intramedullary canal.

Stems allow stability of the spacer and keep it in the joint, and have anti-rotation design. The hinge is manufactured with a dorsal T-shaped design to prevent hyperextension. Stems also have an inclination in relation to the hinge of 15° (PIP) and 30° (MCP), pre-bending angle, to respect the normal position of the fingers and avoid an overload of the silicone material.

The presence of an articulated hinge, between the pivots, allows to restore the finger dorsiflexion movement.

1.1 INTENDED USE

The intended use of the device is the arthroplasty of MCP and PIP joints of the hand. Implants are intended to the treatment of degenerative or inflammatory diseases, dislocations or subluxations for which MCP or PIP joint arthroplasty is necessary, such as:

- Rheumatoid arthritis
- Arthrosis
- Anylocking joints or joints with limited range of motion that did not respond to conservative treatment
- Non-functional joints due to inadequate bone alignment and joint space that cannot be restored by soft tissue reconstruction only
- Joint surfaces destroyed.

The devices represent a permanent implant and are designed to be used as spacers between the bones, to relieve pain and promote greater range of motion in the extremities.

2 CLASSIFICATION

Classification according to Directive 93/42/EEC, Annex IX, rule 8, class IIB.

ATTENTION!
THE DEVICE IS INTENDED FOR USE BY SPECIALIZED ORTHOPEDIC SURGEONS, PROPERLY TRAINED IN PERFORMING HAND ARTHROPLASTY

3.3 USEFUL LIFE OF DEVICE

The device is provided sterile and is intended for single use only. The useful life of the device, therefore, corresponds to its use or to the end of the sterility guarantee, whose date is clearly indicated on the packaging.

Once implanted, the device has a useful life of 5 years.

During this period, it is necessary to carry out periodic checks, regarding the content of section 2.1 INFORMATION TO BE PROVIDED TO THE PATIENT.

1.6 KNOWN CONTRAINDICATIONS AND SIDE EFFECTS

The use of Digitalis device is generally not recommended in the following cases:

- Inadequate muscle-tendon system and skin;
- Inadequacy of the neuro-vascular system;
- Bone demineralization at a significant stage;
- Inadequate bone formation and quantity;
- Ongoing infection;
- Active sepsis;
- The psychological state of the patient is such that the implant is discouraged;
- Patients in childhood or who have not reached skeletal maturity yet.

Use of the device for revision surgeries of implants other than Digitalis is prohibited.

The possibility of developing the following side effects due to the characteristics of the device or possible implant errors should be considered:

- Infection
- Device break
- Pain
- Device displacement

- Hypersensitivity to silicone
- Reaction to foreign body
- Damage to soft tissues, nerves, blood vessels and surrounding tendons
- Loss or reduction of joint mobility
- Synovitis
- Complications of bone cysts.

1.5 SYMBOLOGY

To make the reading of the manual pleasant and clear, the symbols used to manage important warnings for proper and safe use of the device are shown below.

Requirement for proper use

This symbol identifies the presence of information for proper use of the device.

Requirement of information

This symbol identifies the presence of useful and general information, whose reading guides the user towards the conscious use of the device and/or the performance of actions.

Identifies that the product is manufactured, designed and produced in compliance with the safety requirements (RES) of the Medical Devices Directive 93/42/EEC (Class IIB medical device, compliant with classification rule 8, according to Annex IX).

1.6 PRELIMINARY NOTICES

Failure to follow the instructions and precautions below, as well as the rules and precautions described in this user manual, will result in the immediate loss of any warranty on the Digitalis device.

BRM Extremities S.r.l. is not responsible for any damage to people or property as a result of non-compliance with the rules or precautions listed below and generally reported in this user manual.

FAILURE OR NEGLIGENCE IN COMPLYING WITH THE FOLLOWING INSTRUCTIONS MAY CAUSE MALFUNCTION, DAMAGE AND INJURY TO THE USER OR PATIENT

DO NOT USE THE DEVICE UNTIL YOU HAVE FULLY READ AND UNDERSTOOD THIS USER MANUAL

NO DEVICE MODIFICATION IS ALLOWED

USE OF THE DEVICE FOR PURPOSES OTHER THAN THOSE INDICATED IN THIS USER'S MANUAL MAY EXPOSE THE PATIENT AND THE USER TO HAZARD

DO NOT USE THE DEVICE TO REPLACE SIMILAR PROSTHESES FROM ANOTHER MANUFACTURER. THE STEM SIZES MAY NOT BE COMPATIBLE, WHICH COULD LEAD TO MOBILIZATION OF THE DEVICE

THE DEVICE IS DISPOSABLE; RE-USE AND ANY REPROCESSING ACTIVITY IS STRICTLY PROHIBITED.

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Digitalis device is supplied sterile (ethylene oxide sterilization) in a cardboard package. Inside the box there is a double envelope containing the silicone device. The box also contains these instructions for use and traceability labels for medical records.

Digitalis device is sold in 9 different sizes, 4 for the PIP joint and 5 for the MCP joint, as listed below.

REF	Description
DDG3T01001	MCP SPACER DIGITALIS - SIZE 1
DDG3T01002	MCP SPACER DIGITALIS - SIZE 2
DDG3T01003	MCP SPACER DIGITALIS - SIZE 3
DDG3T01004	MCP SPACER DIGITALIS -

BRM Extremities S.r.l.
Sede social - Via Lorenzo Mascheroni 29, 20145 Milão (MI) - Italia
Sede operacional - Via Privata ai Campi 7, 23862 Civate (LC) - Italia
Número IVA - 03416310961
Tel. - 0341 1693087

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