



# BRM BURRS

13 Different cutting edge geometries  
60 mm stem lengths




- Shannon
- Shannon Isham
- Conica
- Wedge
- Cylindrical straight
- Rounded
- Brophy



# BRM BURRS


SHANNON SHORT  $\varnothing$  2.0 x 8 mm

(cylindrical helical teeth)




SHANNON STD  $\varnothing$  2.0 x 13 mm

(cylindrical helical teeth)



SHANNON ISHAM  $\varnothing$  2.0 x 12 mm


(cylindrical straight teeth)



---


CONICAL  $\varnothing$  3.1 x 13 mm

(conical helical teeth)



CONICAL  $\varnothing$  4.1 x 13 mm


(conical helical teeth)



---


CYLINDRICAL STRAIGHT  $\varnothing$  2.0 x 20 mm

(cylindrical straight teeth)



CYLINDRICAL STRAIGHT  $\varnothing$  2.5 x 14 mm

(cylindrical straight teeth)



CYLINDRICAL STRAIGHT  $\varnothing$  3.1 x 22 mm


(cylindrical straight teeth)



# BRM BURRS


**WEDGE Ø 3.1 x 13 mm**

**(helical conical)**




**WEDGE Ø 4.1 x 13 mm**

**(helical conical)**




**CYLINDRICAL ROUNDED Ø 3.1 x 15 mm**

**(cylindrical straight teeth)**




**CYLINDRICAL ROUNDED Ø 4.1 x 15 mm**

**(cylindrical straight teeth)**



**BROPHY Ø 5.0 x 15 mm**

**(helical oval)**



## CODE

## DESCRIPTION

BRB01S62008S	SHANNON SHORT BURR Ø2.0 X 8.0 MM - PACK OF 1 - STERILE
BRB02S62013S	SHANNON STD BURR Ø2.0 X 13 MM - PACK OF 1 - STERILE
BRB03S62012S	SHANNON ISHAM BURR Ø2.0 X 12 MM - PACK OF 1 - STERILE
BRB04S63113S	CONICAL BURR Ø3.1 X 13 MM - PACK OF 1 - STERILE
BRB05S64113S	CONICAL BURR Ø4.1 X 13 MM - PACK OF 1 - STERILE
BRB06S62020S	CYLIN.L STRAIGHT BURR Ø2.0 X 20 MM - PACK OF 1 - STERILE
BRB07S62514S	CYLIN. STRAIGHT BURR Ø2.5 X 14 MM - PACK OF 1 - STERILE
BRB08S63122S	CYLIN. STRAIGHT BURR Ø3.1 X 22 MM- PACK OF 1 - STERILE
BRB09S63113S	WEDGE BURR Ø3.1 X 13 MM - PACK OF 1 - STERILE
BRB10S64113S	WEDGE BURR Ø4.1 X 13 MM - PACK OF 1 - STERILE
BRB11S63115S	ROUNDED BURR Ø3.1 X 15 MM - PACK OF 1 - STERILE
BRB12S64115S	ROUNDED BURR Ø4.1 X 15 MM - PACK OF 1 - STERILE
BRB13S65015S	BROPHY BURR Ø5.0 X 15 MM - PACK OF 1 - STERILE

