

TOOL PRO SCREWS

The choice of compression

Lorem ipsum



Self-drilling and self-tapping
Fully cannulated for guide wire
Anodized titanium alloy

- Fixation of medium and large bone fragments
- Epiphyseal detachments
- Tibio-tarsal arthrodesis
- Ankylosis of ankle and wrist



SHORT THREAD SCREWS *

STERILE	Ø	L
SOS156030CS	5.6	30
SOS156035CS	5.6	35
SOS156040CS	5.6	40
SOS156045CS	5.6	45
SOS156050CS	5.6	50
SOS156055CS	5.6	45
SOS156060CS	5.6	60
SOS156065CS	5.6	65
SOS156070CS	5.6	70
SOS156075CS	5.6	75
SOS156080CS	5.6	80
SOS156085CS	5.6	85

FULL THREAD SCREWS *

STERILE	Ø	L
SOS156035TS	5.6	35
SOS156040TS	5.6	40
SOS156045TS	5.6	45
SOS156050TS	5.6	50
SOS156055TS	5.6	55
SOS156060TS	5.6	60
SOS156065TS	5.6	65
SOS156070TS	5.6	70
SOS156075TS	5.6	75
SOS156080TS	5.6	80
SOS156085TS	5.6	85

STERILE	Ø	L
SOS17040CS	7.0	40
SOS17045CS	7.0	45
SOS17050CS	7.0	50
SOS17055CS	7.0	55
SOS17060CS	7.0	60
SOS17065CS	7.0	65
SOS17070CS	7.0	70
SOS17075CS	7.0	75
SOS17080CS	7.0	80
SOS17085CS	7.0	85
SOS17090CS	7.0	90
SOS17095CS	7.0	95
SOS170100CS	7.0	100
SOS170105CS	7.0	105
SOS170110CS	7.0	110

STERILE	Ø	L
SOS17050TS	7.0	50
SOS17055TS	7.0	55
SOS17060TS	7.0	60
SOS17065TS	7.0	65
SOS17070TS	7.0	70
SOS17075TS	7.0	75
SOS17080TS	7.0	80
SOS17085TS	7.0	85
SOS17090TS	7.0	90
SOS17095TS	7.0	95
SOS17100TS	7.0	100
SOS17105TS	7.0	105
SOS17110TS	7.0	110

* Non-sterile screws available upon request


KIRSCHNER'S WIRE



CODE*	DESCRIPTION	SCREWS
*For the references please contact the manufacturer	K - wire Ø1.6 X 80 mm - pack of 10 - Not Sterile K - wire Ø1.6 X 120 mm with marking at 80 mm - pack of 10 - Not Sterile K - wire Ø1.6 X 150 mm with marking at 80 mm - pack of 10 - Not Sterile	Screws 5.6
*For the references please contact the manufacturer	K - wire Ø2.5 X 80 mm - pack of 10 - Not Sterile K - wire Ø2.5 X 120 mm with marking at 80 mm - pack of 10 - Not Sterile K - wire Ø2.5 X 150 mm with marking at 80 mm - pack of 10 - Not Sterile	Screws 7.0

INSTRUMENTS

CODE	DESCRIPTION
SET.TOOL5670P	Screws instrument box Ø5.6-7.0



BRM EXTREMITIES Srl
Via Lorenzo Mascheroni 29
20145 Milano (MI) Italy
t+39 (02)89415508
www.brm-extremities.com

Headquarter
Via Papa Giovanni XXIII, 9
23862 Civate (LC) Italy
t+39 (0341)1693087
info@brm-extremities.com