



TOOL VIP SCREWS

Verified Intrathread Performance



Self-drilling and self-tapping
Fully cannulated for guide wire
45° chamfered head
Anodized titanium alloy

- Fixation of small and medium-sized bone fragments
- Epiphyseal detachments
- Osteotomy for correction of hallux valgus
- Distal osteotomy of the minor metatarsal for painful plantar callosities
- Ankylosis of ankle and wrist



SCREWS 2.5 *	STERILE	Ø	L	SCREWS 3.0 *	STERILE	Ø	L	SCREWS 4.0 *	STERILE	Ø	L
	SSL125016S	2.5	16		SSL130016S	3.0	16		SSL140016S	4.0	16
	SSL125018S	2.5	18		SSL130018S	3.0	18		SSL140018S	4.0	18
	SSL125020S	2.5	20		SSL130020S	3.0	20		SSL140020S	4.0	20
	SSL125022S	2.5	22		SSL130022S	3.0	22		SSL140022S	4.0	22
	SSL125024S	2.5	24		SSL130024S	3.0	24		SSL140024S	4.0	24
	SSL125026S	2.5	26		SSL130026S	3.0	26		SSL140026S	4.0	26
	SSL125028S	2.5	28		SSL130028S	3.0	28		SSL140028S	4.0	28
	SSL125030S	2.5	30		SSL130030S	3.0	30		SSL140030S	4.0	30
	SSL125032S	2.5	32		SSL130032S	3.0	32		SSL140032S	4.0	32
	SSL125034S	2.5	34		SSL130034S	3.0	34		SSL140034S	4.0	34
	SSL125036S	2.5	36		SSL130036S	3.0	36		SSL140036S	4.0	36
	SSL125038S	2.5	38		SSL130038S	3.0	38		SSL140038S	4.0	38
	SSL125040S	2.5	40		SSL130040S	3.0	40		SSL140040S	4.0	40
					SSL130042S	3.0	42		SSL140042S	4.0	42
					SSL130044S	3.0	44		SSL140044S	4.0	44
					SSL130046S	3.0	46		SSL140046S	4.0	46
					SSL130048S	3.0	48		SSL140048S	4.0	48
					SSL130050S	3.0	50		SSL140050S	4.0	50
					SSL130052S	3.0	52		SSL140052S	4.0	52
					SSL130054S	3.0	54		SSL140054S	4.0	54
					SSL130056S	3.0	56		SSL140056S	4.0	56
					SSL130058S	3.0	58		SSL140058S	4.0	58
					SSL130060S	3.0	60		SSL140060S	4.0	60

FILO DI KIRSCHNER



CODE*

DESCRIPTION

SCREWS

*For the references please
contact the manufacturer

K - wire Ø1.0 X 80 mm - pack of 10 - Not Sterile
K - wire Ø1.0 X 120 mm with marking at 80 mm - pack of 10 - Not Sterile
K - wire Ø1.0 X 150 mm with marking at 80 mm - pack of 10 - Not Sterile

Screws 2.5-3.0

*For the references please
contact the manufacturer

K - wire Ø1.5 X 80 mm - pack of 10 - Not Sterile
K - wire Ø1.5 X 120 mm with marking at 80 mm - pack of 10 - Not Sterile
K - wire Ø1.5 X 150 mm with marking at 80 mm - pack of 10 - Not Sterile

Screws 4.0

INSTRUMENTS

CODE

SET.TOOLVIPP

DESCRIPTION

TOOL VIP Screws instruments box



*Non-sterile screws available upon request