




DIGITALIS



MCP & PIP
JOINTS
ARTHROPLASTY



CE 0426

PRODUCT DESCRIPTION

Elastomeric silicone spacers for metacarpophalangeal (MCP) and proximal interphalangeal (PIP) hand joints arthroplasty.

The Digitalis spacer restores joint mobility and alignment while eliminating pain.

The system is available in different sizes with a pre-flexion angle of 30° for MCP implants and 15° for PIP implants.

ABSOLUTE CONTRAINDICATIONS

- Patient allergy to the implant material
- Acute local or systemic infections
- Severe vascular, neurological, musculo-tendinous, or cutaneous disorders affecting the involved anatomical area
- Insufficient bone structures to ensure adequate implant fixation
- Local bone tumors unless previously treated
- Osteoporosis, osteomalacia, or dysmorphism unless previously treated

CLINICAL INDICATION

- Rheumatoid arthritis
- Osteoarthritis
- Ankylosed joints or joints with a limited range of motion
- Non-functional joints due to inadequate bone alignment and joint space that cannot be restored by soft-tissue reconstruction alone
- Joint surfaces destroyed following traumatic events

- Patients who have not yet reached skeletal maturity
- Severe systemic diseases or metabolic disorders involving the affected anatomical area
- Revision procedures of prosthetic implants manufactured by other companies



For additional information and document downloads

IMPLANTS



PIP	Size
DDG3T02001	1
DDG3T02002	2
DDG3T02003	3
DDG3T02004	4



MCP	Size
DDG3T01001	1
DDG3T01002	2
DDG3T01003	3
DDG3T01004	4
DDG3T01005	5

INSTRUMENTATION

DDGI000001B	Digitalis instruments sterilization box (empty)
DDGI203000	Digitalis Handle
DDGI201003	Multisize Trial MCP
DDGI201004	Multisize Trial PIP
DDGI201001	Digitalis MCP Reamer
DDGI201002	Digitalis PIP Reamer
DDGI202001	Digitalis MCP Rasp
DDGI202002	Digitalis PIP Rasp

DDGI201013	Silicon MCP Trial size 1
DDGI201023	Silicon MCP Trial size 2
DDGI201033	Silicon MCP Trial size 3
DDGI201043	Silicon MCP Trial size 4
DDGI201053	Silicon MCP Trial size 5
DDGI201014	Silicon PIP Trial size 1
DDGI201024	Silicon PIP Trial size 2
DDGI201034	Silicon PIP Trial size 3
DDGI201044	Silicon PIP Trial size 4

KIT DIGITALIS



Distributed by:



BRM EXTREMITIES Srl

Via Lorenzo Mascheroni 29
20145 Milano (MI) Italy
T +39 02 89415508
brm-extremities.com

Headquarter

Via Privata ai Campi 7
23862 Civate (LC) Italy
T +39 0341 1693087
info@brm-extremities.com

D-04A-001-A-EN-REV00