



BRM BURRS

**PERCUTANEOUS
BURRS**

CE 0426




PRODUCT DESCRIPTION

Burrs for percutaneous foot surgery. Helical, cylindrical, conical, or oval burs available in various diameters and lengths for osteotomy and reshaping of the foot bones.

INTENDED USE

The cutting (osteotomy) and/or milling of bone tissue during minimally invasive percutaneous foot surgery procedures.



STAINLESS STEEL

13 DIFFERENT CUTTING EDGE GEOMETRIES

SHANK TYPE 2

Ø 2.35 mm (ISO 1797)

L 60 mm



For additional information and document downloads

- **Shannon Short and Standard**

Cylindrical shape and helical cutting edges

- **Wedge**

Conical shape and helical cutting edges

- **Conical**

Conical shape and helical cutting edges

- **Brophy**

Large, cylindrical cutter with helical cutting edges and rounded tip

- **Shannon Isham**

Cylindrical shape and straight cutting edges


- **Cylindrical Straight**

Cylindrical shape and straight cutting edges

- **Rounded**

Cylindrical straight cutting edges with a rounded tip


DEVICES



BRB01S62008S SHANNON SHORT Ø 2.0 x 8 mm
pack of 1 - sterile



BRB02S62013S SHANNON STD Ø 2.0 x 13 mm
pack of 1 - sterile



BRB03S62012S SHANNON ISHAM Ø 2.0 x 12 mm
pack of 1 - sterile




BRB04S63113S CONICAL Ø 3.1 x 13 mm
pack of 1 - sterile




BRB05S64113S CONICAL Ø 4.1 x 13 mm
pack of 1 - sterile


DEVICES



BRB06S62020S CYLINDRICAL STRAIGHT Ø 2.0 x 20 mm
pack of 1 - sterile




BRB07S62514S CYLINDRICAL STRAIGHT Ø 2.5 x 14 mm
pack of 1 - sterile



BRB08S63122S CYLINDRICAL STRAIGHT Ø 3.1 x 22 mm
pack of 1 - sterile



BRB09S63113S WEDGE Ø 3.1 x 13 mm
pack of 1 - sterile




BRB10S64113S WEDGE Ø 4.1 x 13 mm
pack of 1 - sterile


BRM BURRS

BRM BURRS


DEVICES



BRB11S63115S CYLINDRICAL ROUNDED Ø 3.1 x 15 mm
pack of 1 - sterile



BRB12S64115S CYLINDRICAL ROUNDED Ø 4.1 x 15 mm
pack of 1 - sterile



BRB13S65015S BROPHY Ø 5.0 x 15 mm
pack of 1 - sterile



BRM BURRS

INSTRUMENTATION



STVIPR0001 Beaver handle



STVIPR0002 Straight Elevator



STVIPR0003 Double Rasp



STVIPR0004 Curved Elevator



BRM-MI006-00 Sterilization box for percutaneous instruments

PERCUTANEOUS KIT



Distributed by:



BRM EXTREMITIES Srl

Via Lorenzo Mascheroni 29
20145 Milano (MI) Italy
T +39 02 89415508
brm-extremities.com

Sede operativa

Via Privata ai Campi 7
23862 Civate (LC) Italy
T +39 0341 1693087
info@brm-extremities.com

D-06B-001-A-EN REV0