




Silktoe®



**FIRST METATARSO-
PHALANGEAL (MTP1)
JOINT ARTHROPLASTY**

CE 0426



BRM®
Extremities

PRODUCT DESCRIPTION

Elastomeric silicone spacer for first metatarsophalangeal (MTP1) joint arthroplasty.

The silicone body replaces the joint, restoring the natural dorsiflexion movement of the big toe.

The two metallic rings (grommets) made of titanium alloy act as the interface between the spacer and the bone stumps of the phalanx and metatarsal head.

CLINICAL INDICATION

Moderate or advanced stage hallux rigidus that has not responded adequately to conservative treatments


ABSOLUTE CONTRAINDICATIONS

- Patient allergy to implant materials
- Acute local or systemic infections
- Severe vascular, neurological, musculoskeletal, or skin disorders affecting the anatomical area involved
- Insufficient bone structures to ensure adequate implant fixation
- Local bone tumors if not intentionally treated
- Osteoporosis, osteomalacia, or skeletal deformities if not intentionally treated
- Patients who have not yet reached skeletal maturity
- Pregnant women
- Severe systemic diseases or metabolic disorders involving the anatomical area concerned



For additional information and document downloads





IMPLANTS

REFERENCE	SIZE
BRDNEWP2020L	20 LEFT
BRDNEWP3030L	30 LEFT
BRDNEWP3535L	35 LEFT
BRDNEWP4040L	40 LEFT
BRDNEWP5050L	50 LEFT

IMPLANTS

REFERENCE	SIZE
BRDNEWP2020R	20 RIGHT
BRDNEWP3030R	30 RIGHT
BRDNEWP3535R	35 RIGHT
BRDNEWP4040R	40 RIGHT
BRDNEWP5050R	50 RIGHT

INSTRUMENTATION



BRDNPTHANDLE-QK Silktoe Quick Handle



BRDNPTQK-TS Quick connection type-S

BRDNPTQK-TC Quick connection type-C




BRDNPTCG2030 20-30 Silktoe Cutting Guide

BRDNPTCG3535 35 Silktoe Cutting Guide

BRDNPTCG4050 40-50 Silktoe Cutting Guide

BRDNPTCGMASK Silktoe Cutting Guide holder



BRDNPTI2030P-QK 20-30 Silktoe Metatarsal Grommets Quick Impactor


BRDNPTI3535P-QK 35 Silktoe Metatarsal Grommets Quick Impactor

BRDNPTI4050P-QK 40-50 Silktoe Metatarsal Grommets Quick Impactor

BRDNPTI2030D-QK 20-30 Silktoe Phalangeal Grommets Quick Impactor

BRDNPTI3535D-QK 35 Silktoe Phalangeal Grommets Quick Impactor

BRDNPTI4050D-QK 40-50 Silktoe Phalangeal Grommets Quick Impactor



BRDNPTR2030P-QK 20-30 Silktoe Metatarsal Quick Rasp

BRDNPTR3535P-QK 35 Silktoe Metatarsal Quick Rasp


BRDNPTR4050P-QK 40-50 Silktoe Metatarsal Quick Rasp

BRDNPTR2030D-QK 20-30 Silktoe Phalangeal Quick Rasp

BRDNPTR3535D-QK 35 Silktoe Phalangeal Quick Rasp

BRDNPTR4050D-QK 40-50 Silktoe Phalangeal Quick Rasp

INSTRUMENTATION



BRDNPTTEST20L 20 Left Silktoe Trial Sizer (Green)

BRDNPTTEST30L 30 Left Silktoe Trial Sizer (Orange)

BRDNPTTEST35L 35 Left Silktoe Trial Sizer (Yellow)

BRDNPTTEST40L 40 Left Silktoe Trial Sizer (Blue)

BRDNPTTEST50L 50 Left Silktoe Trial Sizer (Red)

BRDNPTTEST20R 20 Right Silktoe Trial Sizer (Green)

BRDNPTTEST30R 30 Right Silktoe Trial Sizer (Orange)


BRDNPTTEST35R 35 Right Silktoe Trial Sizer (Yellow)

BRDNPTTEST40R 40 Right Silktoe Trial Sizer (Blue)

BRDNPTTEST50R 50 Right Silktoe Trial Sizer (Red)

BRDNEWPTSBOX-QK Silktoe Sterilization Quick Box (empty)

KIT SILKTOE-QK



Distributed by:



BRM EXTREMITIES Srl

Via Lorenzo Mascheroni 29
20145 Milano (MI) Italy
T +39 02 89415508
brm-extremities.com

Headquarter

Via Privata ai Campi 7
23862 Civate (LC) Italy
T +39 0341 1693087
info@brm-extremities.com

D-02A-001-A-EN-REV00