

BRM TOOL PRO SCREWS

THE CHOICE OF
COMPRESSION

CE 0426



PRODUCT DESCRIPTION


“Headless” interfragmentary lag screws for osteosynthesis of large bone fragments. They combine great self-tapping and self-drilling properties with great fracture reduction and compression capabilities, to provide optimal internal fixation.

Available in 2 different diameters and 2 different types of thread (partial and total).

CLINICAL INDICATION

Fixation of medium and large bone fragments:

- epiphyseal detachments
- ankylosis of the ankle
- ankylosis of the wrist
- calcaneal osteotomy



TORX CONNECTION
SELF-DRILLING
AND SELF-TAPPING
TITANIUM ALLOY Ti6AL4V



For additional information and document downloads

Diameter	5.6	7.0
L screw	From 30 to 85 mm	From 40 to 110 mm
Connection	T20	T25
K-Wire	Ø 1.6	Ø 2.5
L thread	Partial thread / Total thread	
Color code	Grey	Sky Blue



STERILE SCREWS 5.6



Short thread	Ø	L
SOS156030CS	5.6	30
SOS156035CS	5.6	35
SOS156040CS	5.6	40
SOS156045CS	5.6	45
SOS156050CS	5.6	50
SOS156055CS	5.6	45
SOS156060CS	5.6	60
SOS156065CS	5.6	65
SOS156070CS	5.6	70
SOS156075CS	5.6	75
SOS156080CS	5.6	80
SOS156085CS	5.6	85



Full thread	Ø	L
SOS156035TS	5.6	35
SOS156040TS	5.6	40
SOS156045TS	5.6	45
SOS156050TS	5.6	50
SOS156055TS	5.6	55
SOS156060TS	5.6	60
SOS156065TS	5.6	65
SOS156070TS	5.6	70
SOS156075TS	5.6	75
SOS156080TS	5.6	80
SOS156085TS	5.6	85

TOOL PRO
SCREWS

STERILE SCREWS 7.0



Short thread	Ø	L
SOS17040CS	7.0	40
SOS17045CS	7.0	45
SOS17050CS	7.0	50
SOS17055CS	7.0	55
SOS17060CS	7.0	60
SOS17065CS	7.0	65
SOS17070CS	7.0	70
SOS17075CS	7.0	75
SOS17080CS	7.0	80
SOS17085CS	7.0	85
SOS17090CS	7.0	90
SOS17095CS	7.0	95
SOS170100CS	7.0	100
SOS170105CS	7.0	105
SOS170110CS	7.0	110



Full thread	Ø	L
SOS17050TS	7.0	50
SOS17055TS	7.0	55
SOS17060TS	7.0	60
SOS17065TS	7.0	65
SOS17070TS	7.0	70
SOS17075TS	7.0	75
SOS17080TS	7.0	80
SOS17085TS	7.0	85
SOS17090TS	7.0	90
SOS17095TS	7.0	95
SOS17100TS	7.0	100
SOS17105TS	7.0	105
SOS17110TS	7.0	110

TOOL PRO
SCREWS

INSTRUMENTATION



UAOI0SA00L	Quick Coupling Handle	for screws Ø 5.6 mm
UAOI0MT000	T-Type Quick Coupling Handle	for screws Ø 7.0 mm
UAOI0MJ100	J-Type Quick Coupling Handle	for screws Ø 7.0 mm



STAIGW016	Guidewire Ø 1.6 mm	for screws Ø 5.6 mm
STAIGD039	Guidedrill Ø 3.9 mm	for screws Ø 5.6 mm
STAIPS056	Tissue Protection Sleeve Ø 5.6 mm	for screws Ø 5.6 mm
STAIGW025	Guidewire Ø 2.5 mm	for screws Ø 7.0 mm
STAIGD046	Guidedrill Ø 4.6 mm	for screws Ø 7.0 mm
STAIPS070	Tissue Protection Sleeve Ø 7.0 mm	for screws Ø 7.0 mm



STAI256150M	Depth Gauge 85 mm	KW 150	for screws Ø 5.6 mm
STAI270000M	Depth Gauge 110 mm	KW 200	for screws Ø 7.0 mm



STAI239GD0P	Cannulated Drill Bit	Ø 3.9 mm	for screws Ø 5.6 mm
STAI246GD0P	Cannulated Drill Bit	Ø 4.6 mm	for screws Ø 7.0 mm



STAI256000K	Basic Countersink	Ø 5.6 mm	for screws Ø 5.6 mm
STAI270000K	Basic Countersink	Ø 7.0 mm	for screws Ø 7.0 mm



STAI100020L	Long Cannulated Screwdriver	BTX20	for screws Ø 5.6 mm
STAI100025L	Long Cannulated Screwdriver	BTX25	for screws Ø 7.0 mm



STAI100020PL	Long Solid Screwdriver	BTX20	for screws Ø 5.6 mm
STAI100025PL	Long Solid Screwdriver	BTX25	for screws Ø 7.0 mm


STOI056070BPL	Metallic Sterilization tray for Instruments Tool Pro
---------------	--

K-WIRES



STOIKW16201503N	K-Wire Ø 1.6 x 200 mm with marker at 150 mm Pack of 3 pcs - Non Sterile	for screws Ø 5.6 mm
STOIKW16201510N	K-Wire Ø 1.6 x 200 mm with marker at 150 mm Pack of 10 pcs - Non Sterile	for screws Ø 5.6 mm
STOIKW25200003N	K-Wire Ø 2.5 x 200 mm Pack of 3 pcs - Non Sterile	for screws Ø 7.0 mm
STOIKW252000010N	K-Wire Ø 2.5 x 200 mm Pack of 10 pcs - Non Sterile	for screws Ø 7.0 mm

KIT TOOL PRO



Distributed by:



BRM EXTREMITIES Srl

Via Lorenzo Mascheroni 29
20145 Milano (MI) Italy
T +39 02 89415508
brm-extremities.com

Headquarter

Via Privata ai Campi 7
23862 Civate (LC) Italy
T +39 0341 1693087
info@brm-extremities.com

D-01A-003-A-EN-REV00