

BRM TOOL SCREWS

MODERATE
COMPRESSION
BEST
FIXATION


CE 0426



PRODUCT DESCRIPTION

“Headless” interfragmentary lag screws for osteosynthesis of small bone fragments. They combine great self-tapping and self-drilling properties with great fracture reduction and compression capabilities, to provide optimal internal fixation.

Available in 4 different diameters and different thread lengths.



CLINICAL INDICATION

- Fixation of small and medium sized bone fragments:
- epiphyseal detachments
 - osteotomy for the correction of hallux valgus
 - distal osteotomy of the second, third or fourth metatarsal for painful plantar calluses
 - ankylosis of the ankle
 - ankylosis of the wrist
 - scaphoid fractures
 - metacarpal and phalangeal fractures

TORX CONNECTION
SELF-DRILLING
AND SELF-TAPPING
TITANIUM ALLOY Ti6Al4V



For additional information and document downloads

Diameter	2.2	2.6	3.2	4.0
L screw	From 10 to 40 mm		From 10 to 60 mm	
Connection	T7		T8	T9
K-Wire	Ø 0.8		Ø 1.2	
L thread	Long thread / Short thread			Long thread
Color code	Blue	Yellow	Green	Fuchsia

IMPLANTS



STERILE SCREWS 2.2



Long thread	Ø	L
STO122010PS	2.2	10
STO122012PS	2.2	12
STO122014PS	2.2	14
STO122016PS	2.2	16
STO122018PS	2.2	18
STO122020PS	2.2	20
STO122022PS	2.2	22
STO122024PS	2.2	24
STO122026PS	2.2	26
STO122028PS	2.2	28
STO122030PS	2.2	30
STO122032PS	2.2	32
STO122034PS	2.2	34
STO122036PS	2.2	36
STO122038PS	2.2	38
STO122040PS	2.2	40



Short thread	Ø	L
STO122011PS	2.2	11
STO122013PS	2.2	13
STO122015PS	2.2	15
STO122017PS	2.2	17
STO122019PS	2.2	19
STO122021PS	2.2	21
STO122023PS	2.2	23
STO122025PS	2.2	25
STO122027PS	2.2	27
STO122029PS	2.2	29
STO122031PS	2.2	31
STO122033PS	2.2	33
STO122035PS	2.2	35
STO122037PS	2.2	37
STO122039PS	2.2	39

STERILE SCREWS 2.6



Long thread	Ø	L
STO126010PS	2.6	10
STO126012PS	2.6	12
STO126014PS	2.6	14
STO126016PS	2.6	16
STO126018PS	2.6	18
STO126020PS	2.6	20
STO126022PS	2.6	22
STO126024PS	2.6	24
STO126026PS	2.6	26
STO126028PS	2.6	28
STO126030PS	2.6	30
STO126032PS	2.6	32
STO126034PS	2.6	34
STO126036PS	2.6	36
STO126038PS	2.6	38
STO126040PS	2.6	40



Short thread	Ø	L
STO126011PS	2.6	11
STO126013PS	2.6	13
STO126015PS	2.6	15
STO126017PS	2.6	17
STO126019PS	2.6	19
STO126021PS	2.6	21
STO126023PS	2.6	23
STO126025PS	2.6	25
STO126027PS	2.6	27
STO126029PS	2.6	29
STO126031PS	2.6	31
STO126033PS	2.6	33
STO126035PS	2.6	35
STO126037PS	2.6	37
STO126039PS	2.6	39

STERILE SCREWS 3.2



Long thread	Ø	L
SOS132010PS	3.2	10
SOS132012PS	3.2	12
SOS132014PS	3.2	14
SOS132016PS	3.2	16
SOS132018PS	3.2	18
SOS132020PS	3.2	20
SOS132022PS	3.2	22
SOS132024PS	3.2	24
SOS132026PS	3.2	26
SOS132028PS	3.2	28
SOS132030PS	3.2	30
SOS132032PS	3.2	32
SOS132034PS	3.2	34
SOS132036PS	3.2	36
SOS132038PS	3.2	38
SOS132040PS	3.2	40
SOS132042PS	3.2	42
SOS132044PS	3.2	44
SOS132046PS	3.2	46
SOS132048PS	3.2	48
SOS132050PS	3.2	50
SOS132052PS	3.2	52
SOS132054PS	3.2	54
SOS132056PS	3.2	56
SOS132058PS	3.2	58
SOS132060PS	3.2	60



Short thread	Ø	L
SOS132011PS	3.2	11
SOS132013PS	3.2	13
SOS132015PS	3.2	15
SOS132017PS	3.2	17
SOS132019PS	3.2	19
SOS132021PS	3.2	21
SOS132023PS	3.2	23
SOS132025PS	3.2	25
SOS132027PS	3.2	27
SOS132029PS	3.2	29
SOS132031PS	3.2	31
SOS132033PS	3.2	33
SOS132035PS	3.2	35
SOS132037PS	3.2	37
SOS132039PS	3.2	39
SOS132041PS	3.2	41
SOS132043PS	3.2	43
SOS132045PS	3.2	45
SOS132047PS	3.2	47
SOS132049PS	3.2	49
SOS132051PS	3.2	51
SOS132053PS	3.2	53
SOS132055PS	3.2	55
SOS132057PS	3.2	57
SOS132059PS	3.2	59

STERILE SCREWS 4.0



Long thread	Ø	L
SOS140010PS	4.0	10
SOS140012PS	4.0	12
SOS140014PS	4.0	14
SOS140016PS	4.0	16
SOS140018PS	4.0	18
SOS140020PS	4.0	20
SOS140022PS	4.0	22
SOS140024PS	4.0	24
SOS140026PS	4.0	26
SOS140028PS	4.0	28
SOS140030PS	4.0	30
SOS140032PS	4.0	32
SOS140034PS	4.0	34



Long thread	Ø	L
SOS140036PS	4.0	36
SOS140038PS	4.0	38
SOS140040PS	4.0	40
SOS140042PS	4.0	42
SOS140044PS	4.0	44
SOS140046PS	4.0	46
SOS140048PS	4.0	48
SOS140050PS	4.0	50
SOS140052PS	4.0	52
SOS140054PS	4.0	54
SOS140056PS	4.0	56
SOS140058PS	4.0	58
SOS140060PS	4.0	60

INSTRUMENTATION



UAOI0MA000 Quick Coupling Handle



UAOI0SA00S AO Quick Coupling Silicone Handle



STOIGW012 Wire guide Ø 1.2 mm



STOI232000M Depth gauge 60 mm



STOI217008P Cannulated drill bit Ø 1.7 mm KW 0.8 for screws Ø 2.2 - 2.6 mm

STOI222000P Cannulated drill bit Ø 2.2 mm KW 1.2 for screws Ø 3.2 mm

STOI227012P Cannulated drill bit Ø 2.7 mm KW 1.2 for screws Ø 4.0 mm



STOI22200BK Basic countersink Ø 2.2 - 2.6 mm for screws Ø 2.2 - 2.6 mm

STOI23200BK Basic countersink Ø 3.2 mm for screws Ø 3.2 mm

STOI24000BK Basic countersink Ø 4.0 mm for screws Ø 4.0 mm



STOI108007CS Slim cannulated screwdriver BTX7 KW 0.8 for screws Ø 2.2 - 2.6 mm

STOI100008CS Slim cannulated screwdriver BTX8 KW 1.2 for screws Ø 3.2 mm

STOI112009CS Slim cannulated screwdriver BTX9 KW 1.2 for screws Ø 4.0 mm



STOI100007P Solid screwdriver BTX7 for screws Ø 2.2 - 2.6 mm

STOI100008P Solid screwdriver BTX8 for screws Ø 3.2 mm

STOI100009P Solid screwdriver BTX9 for screws Ø 4.0 mm



UKWI00120T Tube for K-wire sterilization

STOI022040BPH Sterilization Tray

K-WIRES



STOIKW08120803N K-Wire Ø 0.8 x 120 mm with marker at 80 mm for screws Ø 2.2 - 2.6 mm
Pack of 3 pcs - Non Sterile

STOIKW08120810N K-Wire Ø 0.8 x 120 mm with marker at 80 mm for screws Ø 2.2 - 2.6 mm
Pack of 10 pcs - Non Sterile

STOIKW12120803N K-Wire Ø 1.2 x 120 mm with marker at 80 mm for screws Ø 3.2 - 4.0 mm
Pack of 3 pcs - Non Sterile

STOIKW12120810N K-Wire Ø 1.2 x 120 mm with marker at 80 mm for screws Ø 3.2 - 4.0 mm
Pack of 10 pcs - Non Sterile

KIT TOOL



Distributed by:



BRM EXTREMITIES Srl

Via Lorenzo Mascheroni 29
20145 Milano (MI) Italy
T +39 02 89415508
brm-extremities.com

Headquarter

Via Privata ai Campi 7
23862 Civate (LC) Italy
T +39 0341 1693087
info@brm-extremities.com

D-01A-001-A-EN-REV00